HELMUT ZEPF

www.zepf-dental.com

REF 85.285.10 - 85.290.10 / 85.327.00 - 85.327.30

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
www.zepf-dental.com/anwenderinformationen

**User Information
on HELMUT ZEPF products**

www.zepf-dental.com/en/user-information

Trade Name: **Sterile Container Systems**
REF 85.285.10 - 85.290.10 / 85.327.00 - 85.327.30

Dear Customer!

With the purchase of this sterilization container you have chosen a high quality product, the correct handling and use of which are described in the following.

In order to keep risks and unnecessary burdening of the patient, the user and third parties as low as possible, we request that you carefully look through these instructions for use and keep them for further reference.

Please read carefully the warning notices identified by this symbol. Improper use of the products can lead to serious injuries to the patients, the users or third parties.

Notice

- The processing of medical devices must be complied with the national regulations and standards.
- In case of patients who have or are suspected of having Creutzfeldt-Jakob disease (CJD) or any kind of variations of this disease, the currently relevant national regulations must be applied concerning the processing.

Warning notices and precautionary measures

- Only container tubs and lids in perfect condition with an undamaged seal and with an undamaged filter system or inserted filter may be used for sterilization! Only one filter per filter holder may be used.

- Disposable paper filters are preferable to reusable filters.
- A washing and cleaning by a machine is preferable as it performs a more effective result. Using a machine obtains a greater safety by the cleaning and disinfection process.
- Alkaline cleaning agents (pH >10) are not suitable for all materials. The Robert Koch Institute point out some potential problems caused by increased abrasion by aluminium, silicone elastomeric, adhesive connections, solder connections of silver and stannic, sealing materials, plastic coatings, glass fibre light guides and visual surfaces with antireflection features.
- Never use metal brushes, metal scouring pads or abrasive cleaning agents for manual cleaning.
- Only clean and low microbiological contaminated sterilization containers allow a successful sterilization.
- Defective products must run through the entire sterilization reprocessing before returning them for repair or reclamation. The proof of decontamination must be enclosed with the return shipment.

The form is available for download on our website „https://www.zepf-dental.com/IFU/080300008_GB.pdf“.

- The lid should not be loaded without the bottom or used as storage area, to avoid deformation of the lid.
- If sterilization containers are deformed in any way during the sterilization process, there is a fundamental risk of non-sterility. In this case, the entire batch must not be used and the cause must be determined immediately (check the sterilization protocol; check the sterilizer and the other sterilization containers; search for the cause by performing functional tests on the damaged sterilization container).

- If the user has a sensitivity to silicone or aluminum, unprotected contact can cause an allergic reaction.
- Not to be used for implants.
- Changes to the product in the phases downstream of production are not planned.
- Never use cleaning solutions with bleach (such as sodium hypochlorite), because these may cause intensive corrosion.

RESTRICTION ON REPROCESSING

- Repeated or frequent processing as per this instruction has only minor effects on the containers' service life.
- The service life of a sterilization container is mainly determined by wear and damage through usage.
- In case of a proper use on an average of four times a week, arises a service life of containers and baskets of up to 10 years.

1 AREA OF USE

Helmut Zepf sterile container systems combine approved filter technologies, tested materials and design properties into a reliable container system. These are reusable container systems which offer a wide selection on dimensions and equipment in order to provide an effective packaging, storage and transportation of the surgical instruments for sterilizing. The container systems are optimally suitable for fractionated vacuum procedure.

1.1 INTENDED PURPOSE

The Helmut Zepf sterile container systems are intended to be loaded with medical devices for sterilization. The sterilization and storage of the encased products are enabled and insured until the next use.

The containers are available depending on models with perforated and non-perforated bottoms and perforated lids. The 1/1, 3/4 and 1/2 containers are also available with safety lids / covers.

1.2 COMBINATION PRODUCTS

Sterile container systems are a combination of sterile containers, wire / screen baskets and filters. Accessories can also be used for the container systems. For each container size should be used a suitable sized basket. The possible combinations of various container designs are shown in the following.

Standard container

The lids and if applicable, the bottoms are equipped with filter holders below / above the perforations. Before sterilization, disposable paper filters or permanent filters must be inserted into the filter holders.

If necessary, a safety cover can also be placed on the lid of the 1/1, 1/2 and 3/4 sized standard containers. This protects the sterile container from contamination during storage or transport. The new standard model of sterile container has this safety cover already integrated. The safety lid may only be fitted for transportation and must not be used for reprocessing.

Baskets

There are suitable wire baskets in various heights with compatible lids and screen baskets with compatible lids and feet for each container size.

Security seals

Security seals are attached to the closures outside by guiding the seal through the opening of the spring lock system and joining the seal. The seal will break by opening / flipping the latches.

Trade Name: **Sterile Container Systems**

Information for the use of paper filters

- Paper filters are intended for one use only.
- Paper filters are produced according to ISO 11607-1.
- No adhesives may be applied to the paper filters (e.g. to document cycles), since adhesives may contain harmful substances.

Furthermore, adhesion destroys the germ barrier.

The paper filters must have the correct measurements to cover the perforation completely in the container lid / bottom.

Information for the use of permanent filters

- PTFE filters are intended for multiple uses (up to 1200 cycles).
- No adhesives may be applied to the permanent filters (e.g. to document cycles), since adhesives may contain harmful substances. Furthermore, adhesion destroys the germ barrier.
- In case the filter gets particularly contaminated than it must be removed and cleaned.

The permanent filters must have the correct measurements to cover the perforation completely in the container lid / bottom.

2 HANDLING

2.1 GENERAL

The Helmut Zepf sterilization containers are made of an aluminium alloy and its surface is anodized to protect against corrosion. Aggressive cleaning products, metal brushes or scouring pads could permanently damage this surface and therefore may not be used. If this instruction is not followed than it leads to the exclusion of the warranty.

The sterilization containers may only be handled by trained / instructed personnel in order to prevent damage to the containers, closures, gaskets and sterile filters / cassettes. Direct patient contact is not intended.

Sterilization indicators and colored identification labels provide information about the content, place of use and state.

Suitable measures (such as security seals or process indicators) must be taken according to normative requirements and recommendations to ensure that sterilized and unsterilized sterilization containers cannot be mistaken for one another.

Only intact sealing ensures that the sterilization container has not been opened without permission.

2.2 PREPARATION FOR CLEANING

1. Separate the container bottom and lid
2. Remove the contents of the container (basket, instruments, etc.)
3. Remove the filter holders / cassettes from the inside of the lid and if applicable from the bottom part (in the case of containers with perforated bottom)
4. In case of disposable paper filter: dispose the disposable filter
5. Remove the disposable seals and indicator labels

Note: All paper filters are disposable filters and must be replaced after each use of the container.

Note: The lid should not be loaded without the bottom, e.g. with the bottom or used as storage area to avoid deformation of the lid.

2.3 IMPLEMENTING OF A NEW / UNUSED CONTAINER

The container must be thoroughly cleaned before it is used for the first time.

- The container must be pre-processed in a validated cleaning / washing machine and disinfection process.
- A neutral cleaning agent should be used for this purpose in the machine.
- Once the pre-processing in the cleaning / washing machine and disinfection is completed, the products must be steam-sterilized at 134° C in a fractionated steam sterilization process.

- Furthermore, all movable parts on the container must be treated regularly with approved instrument maintenance oil.
- Suitable new filters will have to be inserted after the cleaning process (see 2.6 Filter change).

2.4 CLEANING AND DISINFECTION

Note: Improper washing and disinfection may lead to corrosion and stress fractures. For that reason, the specifications of the washing / disinfection product manufacturer must be considered. The cleaning agent must be free of sodium, alkaline and carbonate and it must have a neutral pH value and / or approved for the treatment of anodized aluminium by the chemical's manufacturer. The water used for this purpose must be in conformity at least with the quality recommendation of the manufacturer of the washer-disinfector for the proper operation of the equipment. The following essential specifications must be followed during processing the containers:

The container must be washed and disinfected before the first use and after each further use.

2.4.1 MECHANICAL CLEANING

Contamination which cannot be removed in the usual washing procedure, (such as adhesive labels, indicator strips and markings) can be removed with an eloxal-cleaning agent. After this special treatment, the products must be processed as described in the following:

- Neutral or other suitable cleaning and disinfection solution are to be used which are explicitly approved for the treatment of aluminium products. If necessary, the products must be checked for the suitability for the relevant procedure.
- Use only neutral cleaning agent for aluminium containers, which are specifically approved for the washing of aluminium containers by the manufacturers. For the proper dosing please refer also to the manufacturer's specifications. These products are also suitable for cleaning surgical instruments by optimising the program.
- If using neutralisation agent then the products must be checked for the suitability for aluminium.
- Low-salt water should necessarily be used for the final rinsing.
- The cleaning devices and inserts must be suitable for the processing of containers and lids. This applies particularly for the correct positioning in the loading inserts to ensure an adequate and unobstructed rinsing, draining of the applied media flow and the drying of all container and lids.
- Containers and lids may not be cleaned and disinfected in a closed / assembled state.
- Attention must be paid when loading the machine to ensure a sufficient media flow during the process.
- The container bottom must be placed into the washing machine with the opening downward in order to prevent the accumulation of water and to ensure that the used media flows off adequately.
- The container lid must be washed with the inside facing downward and the latches / closures folded outward.
- The containers and their accessories, without any visible residues, will be removed from the media after completing the mechanical cleaning and disinfection process.

Should there be still any residues detected then the position of the containers and accessories in the device should be rechecked and possibly changed.

The following cleaning procedure and adjustment has been validated in a suitable cleaning and disinfection device (Miele PG 8528) by us and can be used for the cleaning:

- 1 min. pre-cleaning with cold (< 40 °C) water (drinking water quality according to EC Directive 98/83/EC)
- 3 min. cleaning with Mediclean 0.5% (Dr. Weigert) at 45° C
- Neutralization with demineralized water (DI-water)

Trade Name: MD Sterile Container Systems

2.4.2 MANUAL CLEANING / DISINFECTION

- As far as possible mild, neutral cleaning agents should be used for the aluminium container and lids or chemical products which are specifically approved also for the treatment of aluminium by the manufacturer. If necessary, the products must be checked for the suitability for the relevant procedure. A suitable soft sponge should be used for manual cleaning. Do not use scouring pads because they destroy the surfaces and thereby the passivation and that leads to the loss of the warranty rights.
- After washing, a thoroughly rinsing with suitable low-salt water (such as demineralized water) and sufficient drying is necessary.
- Do not use any metal brushes or abrasive cleanser.
- In the case of the PTFE filters, manual cleaning is only carried out in the event of heavy contamination of the filter. The filter is removed from the container and carefully cleaned. Only hospital-approved cleaning agents for containers and surgical instruments are used.
- Finally, disinfection is to accomplish according to the respective hygiene requirements.

2.5 INSPECTION, MAINTENANCE AND TESTING

The sterilization containers must be inspected for their functionality before each use. Damages on the closures, gasket, filter holders, and filters as well as bent or dented parts indicate for the need of repair of the sterilization container and may not be used.

Do not use defective sterilization containers.

- The durability of the gasket is up to 500 sterilization cycles. After that, the gaskets will have to be checked out with particular care and have to be changed if needed. The operator is responsible that the amount of the conducted cycles of a product can be checked.
- The sealing has to be checked before every use.
- All movable parts on the container must be treated with approved instrument maintenance oil.
- If there are any damages detected on the gaskets, then they have to be replaced immediately.
- The gaskets should not be treated with spray, oil or solvents. It is enough for cleaning and maintenance to wipe occasionally with a moist cloth.
- If there are any damages detected on the sterilization containers then they have to be inspected, repaired or replaced if necessary.
- The sterilization containers may be maintained and repaired only by qualified persons. Do not try to repair yourself the gaskets or attachments in order not to compromise the safe use of the container, unless it applies to containers which have a special replaceable silicone gasket.
- The sterilization containers may be returned for maintenance or repair to Helmut Zepf Medizintechnik GmbH or to an authorized repair service.
- We recommend the use of our medical white oil **24.950.10** for the care of all moving parts on the container. Also available in the practical oil pen **24.950.12**.
- The following replacement parts are available from Helmut Zepf Medizintechnik GmbH:

- Paper filter with indicator 215 x 145 mm (100's.) for 1/1-, 3/4-, 1/2- and mini containers **85.320.01**
- PTFE (Teflon) permanent filter 215 x 145 mm (2's) **85.320.02**
- Paper filter with indicator 235 x 118 mm (100's.) for 1/1-, 3/4-, 1/2- and mini containers **85.327.50**
- PTFE (Teflon) permanent filter 215 x 95 mm (2's) **85.327.60**
- Colored Identification labels **85.312.00**
- brown **85.312.01**
- blue **85.312.02**
- yellow **85.312.03**
- green **85.312.04**
- purple **85.312.05**
- orange **85.312.05**

- red **85.312.07**
- black **85.312.08**
- silver **85.312.09**
- Plastic Security Seals (100's) **85.327.70**

2.6 CHANGING OF FILTER

After changing the filter, the filter holder has to be placed by pressing into its correct position with an audible snap.

Helmut Zepf lids may only be used with Helmut Zepf filter holders.

- Disposable paper sterile filters must be reinserted before every re-sterilization.
- Suitability and proper fit are only guaranteed if Helmut Zepf filters are used.
- Warranty services can only be accepted if original Helmut Zepf filters are exclusively used.
- PTFE filters have been tested for usage duration of 1200 cycles and must be replaced afterwards.

Attention

Combine only original Helmut Zepf component parts such as lids, bottoms, filters, gaskets, cassettes and filter holders with each other in order to avoid putting at risk the leak tightness and germ barrier. Otherwise, Helmut Zepf does not undertake any warranty.

3 STERILIZATION

Helmut Zepf sterilization container systems have been validated with the following sterilization parameters:

- Method: 3 x pre-vacuum steam sterilization
- Temperature: 134° C (273° F)
- Half cycle: 2.5 minutes
- Holding time: 5 minutes
- Drying time: 20 minutes
- Loading: Standard medical instruments (scissors, clamps, forceps) and textiles

The validation covers the following containers of the Helmut Zepf sterilization container system:

- Full-size container, size 1/1
- 3/4-size container
- 1/2-size container
- Dental container
- Mini container

The validation applies to all sterilization containers with the following references: **85.285.10 - 85.290.10 / 85.327.00 - 85.327.30**

Important: We recommend a mechanical cleaning to achieve a more effective result.

Attention

ZEPF sterilization containers have been tested and validated for steam sterilization with moist heat in a fractionated vacuum procedure.

3.1 LOADING OF THE CONTAINER

The total weight of the load of a container should not exceed the following loading size according to DIN EN 868-8. Otherwise a satisfying sterilization result cannot be ensured.

Modell	Max. loading in kg
1/1 (Full-) Size Container	10
3/4 Size Container	7
1/2 Size Container	5
Dental Container	1.8
Mini Container	1

Trade Name: **Sterile Container Systems**

In accordance with national legislation the loading limit can differ from the details listed above.

In case of loading with textiles, please pay attention that the pieces of laundry or folded textiles are in vertical position. It should still be easily possible to slide an open hand between the pieces of laundry by a fully loaded container. The sterilization containers were validated for lumen instruments with diameter from 1.2 mm on to a length of 23.5 mm.

Attention

The sterilization of various container loadings and the configurations must be determined by the responsible personnel.

Endoscopes, instruments with lumens, compressed air or mains-powered units and instruments with cannulas must be prepared for the sterilization according to the specifications of the manufacturer. Small baskets, trays or other accessories, especially the ones with lids or flaps, should only be used together with sterilization container systems, if these were specifically designed and tested for this purpose.

Load limit: For standard containers a loading limit of 20 mm from the top edge of the bottom must be taken into account. For dental and mini containers a loading limit of 5 mm from the top edge of the bottom applies.

Please refer to catalog for further information about Helmut Zepf products, baskets and trays.

Attention

Using water-resistant inserts (such as plastic) may cause remaining condensate inside of the container. Instead, use other suitable mats or holders. Check the integrity of the inserted filter and the proper sealing of the filter holder. Always use the locking mechanism to attach the container lid to the bottom. Ensure that the closures are sealed properly before taking the container for sterilization.

3.2 POSITION IN STERILIZER

The containers are designed so that they can be used in any commercially available large sterilizer for the sterilization with moist heat. Keep in mind that heavy containers are to be positioned at the bottom of the sterilization chamber.

The containers can be stacked easily and safely on top of each other due to their design, without slipping during the sterilization procedure. Stacking is only recommended for sterilization cycles operating with a fractionated vacuum system. The maximum stacking height should not exceed 46 cm in order to ensure an effective air removal and steam penetration. The instructions of the manufacturer of the sterilizer must be followed.

Attention

Pay attention to the followings during sterilization: Never wrap the container in an additional outer packaging. Never cover the perforation are on the lid and bottom with foil packaging or something similar because this prevents the air and steam flow in the container. As a result, there would be a container deformation caused by vacuum due to insufficient pressure compensation and the sterility of the loading cannot be guaranteed. Disposal containers may not be sterilized in closed condition. Always carry the sterile container by the carrying handles and never by the lid during loading and / or unloading the sterilizer as well as during transport.

3.3 SEQUENCE CONTROL

- Operate the loaded sterilizer for the selected sterilizer cycle according to the specifications of the sterilizer manufacturer (referring to temperature and sterilization time). The validation results are in the process to be considered.
- The container should cool down completely on the sterilization cart to avoid condensation in the container.
- The sterile goods must be evaluated and approved after each sterilization according to internal directives and validation results. This is consequently conducted by employees with special knowledge, qualification level 1.

4 STORAGE, TRANSPORT AND DISPOSAL

4.1 STORAGE

Please refer to DIN 58953-9 (Use of sterilization containers) for the storage period of medical devices in sterilization containers. The storage time generally depends on the storage conditions and must be determined by the responsible hygiene specialists. In case of extremely high requirements for asepsis or deviations on the specified storage conditions are shorter storage periods or additional packaging to be used.

Storage conditions:

- Temperature: 15 - 26° C
- Humidity: 30 - 50%
- Air pressure: 500 - 1060 hPa (normal atmosphere pressure)

Different container loadings, storage periods and storage conditions are the responsibility of the hygiene specialists.

Helmut Zepf sterile containers have been tested for a storage period of 12 months by applying of *Bacillus subtilis*, *Aspergillus brasiliensis* and *Candida albicans*. Therefore, depending on the storage conditions a storage period of up to 12 months can be assured. Thereby the containers have to be stored under protected circumstances (e.g. in closed cabinets), protected from dust, clearly, dryly and free of vermin.

4.2 TRANSPORT

The sterile containers should be transported only by the intended carrying handles.

Filters must be protected from damage of any kind during transportation, especially perforations.

4.3 DISPOSAL

For disposal, the products must be free of potentially contaminated material. For this purpose, the product may have to undergo reprocessing.

If there are sharp edges, they must be disposed of in such a way as to avoid endangering people.

5 MATERIALS

The sterilisation containers are made of anodized aluminium alloy and the accessories are made of stainless (surgical) instrument steel.

6 STANDARDS APPLIED

The following standards were considered to ensure the safety of the sterilization containers during handling:

Standard	Description of the standard
DIN EN 868-2	Packaging for terminally sterilized medical devices – Part 2: Sterilization wrap – Requirements and test methods
DIN EN 868-8	Packaging for terminally sterilized medical devices – Part 8: Re-usable sterilization containers for steam sterilizers conforming to EN 285 – Requirements and test methods

Trade Name: **MD Sterile Container Systems**

Standard	Description of the standard
DIN EN ISO 11607-1	Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems
DIN 58952-2	Sterilization – Transport baskets for sterile barrier systems – Part 2: Sterilizing baskets made of metal
DIN 58952-3	Sterilization – Transport baskets for sterile barrier systems – Part 3: Instrument trays for sterilizing goods made of metal
DIN 58953-9	Sterilization; Sterile supply – Part 9: Use of sterilization containers
DIN EN ISO 14937	Sterilization of health care products – General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices
DIN EN ISO 17665	Sterilization of health care products – Moist heat – Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices
DIN EN 285	Sterilization – Steam sterilizers – Large sterilizers

An independent and accredited testing laboratory performed tests in order to ensure the reliability of the sterility. The purpose of these tests was the validation of a sterilization procedure using steam / moist heat by reusable Helmut Zepf sterile container systems. Based on the results, we indicate therefore the specified sterilization procedure on page 1 of this user manual.

10 WARRANTY

These sterilization containers have been manufactured from high quality materials and have been subjected to quality control checks before release to the market. Nevertheless, revert to the address mentioned below in case any errors should arise.

THE HELMUT ZEPF COMPANY ACCEPTS NO LIABILITY IN DEMONSTRABLE CASES OF VIOLATION OF THIS CUSTOMER INFORMATION.

RECOMMENDED CARE PRODUCTS

Medical white oil in the 50 ml dosage bottle and in the oil pen, 12 ml

Physiologically harmless according to the German pharmacopoeia (DAB 10).

24.950.10

24.950.12

RELEASE VERSION – SEE FOOTER

Explanation of symbols			
	Consult instructions for use	UDI	Unique device identifier
	Manufacturer	QTY X	Quantity
	Date of manufacture	HIBC	HIBC Barcode
REF	Catalog number	CE	CE marking
LOT	LOT no. (batch no.)		Caution
	Non-sterile	MD	Medical device