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Trade Name: MD Permanent Filter REF 85.327.60 for MD Sterilization Containers REF 85.327.00 / 85.327.10 / 85.327.15 / 85.327.20 / 85.327.30

Dear Customer!

By purchasing the permanent filters art. no. **85.327.60** you have purchased a high quality product that is described below.

GENERAL

The new permanent filters are a further step in the direction of more safety and economic efficiency by the preparation respectively the sterilization of medical products.

The permanent filters offer the following advantages over disposable filters and previous long-term filters made out of textile material:

- Number of sterilization cycles > 1200
- Suitable for all sterilization container models.
- High economic efficiency through reduction of consumer costs.
- Time saving due to easy filter control and prevention of filter changes.

USAGE / HANDLING

- Handling, respectively insertion of the permanent filter is analogous to disposable filters.
 1. For this purpose, open the internal filter holder in the lid and in the bottom.
 2. If present, remove the filter from there.
 3. Insert a permanent filter.
 4. Fix the permanent filter with the filter holder.
- The permanent filters are only to be handled from instructed, respectively trained persons in order to prevent damage.
- The permanent filter is to be checked before each and every sterilization cycle for damage.
- After the permanent filter is inserted into the filter holder, the filter holder must be depressed so that it can be heard that it clicks into place.

CLEANING

- A mechanical cleaning can be effected without removal of the filter. For cleaning with removal of the filter, please take care that the filter is washed separately. Furthermore, the perforation in the middle of the round filters may not be damaged.
- A manual cleaning has to be effected only if the filter is strongly contaminated. Remove the filter from the container and clean it carefully. Hereby, only detergents approved by the hospital are to be used for the cleaning of containers and surgical instruments.

MATERIAL

PTFE (Teflon®)

The process for filtration of various media with PTFE filters is known for a long time and is now made possible for usage in sterilization containers.

TESTING / EXAMINATIONS

Certification for compatibility of permanent filters was carried out by an independent and accredited laboratory.

Their results are as follows:

The PTFE permanent filters used are suitable to prevent microbial growth according to the requirements of DIN 58953 Part 6, Chapter 2.14.

You may receive the respective report from us upon request.

SPECIAL INFORMATION

Close attention should be paid to the following points in order to guarantee flawless functioning of the permanent filters:

- Never puncture the sterile filters!
- Prevent the placement of sharp and pointed objects on top of the sterilization container.

Prevent the invasion of bacillus through:

- Spilling of dirt, cleaning, or disinfection solutions
- Penetration of contaminated external condensate (i.e. through the spaces by stacking of hot containers)
- Rubbing in of dust and dirt in the filter when wiping off the container
- Long term storage

DISPOSAL

After usage, the permanent filters can be disposed of in the same manner as hospitals do. No special measures are necessary.

GUARANTEE

The permanent filters are manufactured from a high quality material and undergo a strict quality control before delivery. In the event that discrepancies should occur, please contact us directly.

Explanation of symbols	
MD	Medical device
REF	Catalog number
CE	CE marking