ZEPF Kronenspreizzange

ZEPF Crown-Spreading Pliers

REF 19.277.01Z

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Dear valued Customer!

We thank you very much for having chosen this high quality instrument from Helmut Zepf Medizintechnik GmbH.

To ensure that you are satisfied with the instrument for a long period, please find a few instructions for use below.

- 1) First separate the crown completely using the crown-spreading pliers. Please ensure that the crown has been split deep enough. Sharply defined edges with the separator (as cylindrical as possible) give you secure adaption in the slot and greatly facilitate subsequent spreading.
- 2) Now the tip of the crown-spreading pliers is inserted into the slot from occlusal, buccal or lingual, depending on the access option. The crown is getting spreaded by compressing the plier.
- 3) The good anatomical curvature of our crown-spreading pliers allows their use in all quadrants and thus guarantees safe and correct use.
- 4) Careful, directed and slow spreading of the crown now allows widening of the separating gap. This should be at least 1 mm. Compared to normal crown removers, do not apply leverage! This way no pressure is applied to the root, the adjacent tooth or the jaw. The point of rotation is at the tip of the instrument and protects the root underneath.
- 5) Then the crown can be removed without difficulty.

An Allen wrench is provided to position the tips in parallel. If the screw becomes loose during cleaning, the working end can be repositioned after cleaning and before sterilization.

Please note that the wrench is not intended for autoclaving.

CLEANING AND REPROCESSING

For all cleaning processes, please consider our application guidance of autoclavable instruments and accessories according to DIN EN ISO 17664.

See document “0703091007_REPROCESSING_HZ_BA_e“, downloadable from our website www.zepf-dental.com.