
# Benex<sup>®</sup> II

www.benex-dent.com

**ZEPF REF 12.300.08**

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 **Dear Customer!**

We thank you very much for having chosen this high quality instrument from Helmut Zepf Medizintechnik GmbH. To get the best from your instrument, please follow the instructions below:

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 After efficacious anesthesia: Cutting periodontal fibres (Sharpey fibres) in the sulcus by using the **ZEPF** Periostome **Xiod**.

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 Strong, large roots must be loosened / luxated by axial movements within 30 sec., using a slim elevator / twister (**Xiod** from **ZEPF**). Without using transversal movements. In case of multi-rooted teeth, the roots are divided and extracted separately.

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 Drilling with the diamond coated twist drill (REF 12.300.30 / REF 12.300.35) should be in the axis and center of the root fragment. It should be approx. 7 mm in the hard tissue, deeper drilling will not be necessary. Drilling is performed with sufficient cooling, preferably with a physiological, sterile saline solution. In order to remove drilling chips more easily, the drill should be used intermittently, i.e. with inward and outward movements.

**Recommended rotation speed (rpm.): 500-700 | REF 12.300.30 max. 3000 | REF 12.300.35 max. 2200**

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
 Depending on the case situation, the short or long extraction screw with the corresponding diameter is then selected and carefully inserted into the tooth root by using the screwing support.

**Note:** Care must be taken to ensure that the screwing support is pressed into the extraction screw. This is the only way to ensure a sufficient clamping effect for safe transport into the oral cavity.

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The Extractor is positioned on the adjacent crowns: The opening of the round, revolvable segment plate is adjusted in vestibular direction ensuring a good view of the extraction screw. After the traction rope has been hooked into the extraction screw, it is guided over the reverse roller and fixed to the hook of the fully tightened extraction slide.

Under slight traction – so that the rope does not hang out – the instrument is placed on the adjacent teeth by turning the hand screw. When positioning the segment plate, make sure that the screw and the traction rope have the same axial direction.

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
Once the extractor is positioned properly, the extraction is carried out by turning the hand screw. In case of strong, long roots the periodontal fibres have to be pre-stretched during 30 seconds by applying a sub-maximum traction. The extractor is positioned in such a way that it supports itself under tension. This position should be maintained until the periodontal fibers give way after about 30-40 seconds. Subsequently, the periodontal fibers are prestressed again by turning the hand screw. This procedure is repeated until the tooth root is completely extracted.

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**The Free-End Situation:** In order to guarantee an optimized, vertical power-transmission, a right-angle level has to be available in order to apply the Benex® Extractor. With the quadrant support [REF 12.300.80](#) the supporting level can be corrected and/or a bad level can be amended. This is often necessary with Benex® extractions in the mandible and with palatal roots in the maxilla. **Caution:** The use of the quadrant support is also indicated if a neighboring tooth of the tooth to be extracted is crowned. In this case, the quadrant support serves as an even more optimal pressure distribution for the crown.

- A:** The quadrant support is filled with a fast-setting silicone or alginate. In the area of the buccal opening, a recess is modeled.
- B:** Place the quadrant support with opening in the root extraction area, at right angles with the root axis.
- C:** The quadrant support is removed when the silicone or alginate begins to set. The complete setting happens extra-orally.
- D:** Root luxation: Afterwards insertion of the Benex® Screw and repositioning of the quadrant support and the Benex® Extractor.
- E:** Extraction according to steps 5 and 6.

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

In order to achieve an optimal cleaning result, the Benex® Extractor has to be completely disassembled before subjecting it to the cleaning process. The disassembly can be carried out as follows: First of all, remove the "knurled screw" from the extractor by turning the knurled screw anti-clockwise off the thread.

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
In a second step the "swivel segment" is removed from the base plate by turning it anti-clockwise or clockwise. Then remove the "POM segment plate" from the "swivel segment" by pulling it slightly.

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
Step 3 is to remove the "cross pins" from the base plate in a slight pulling or pushing movement. **IMPORTANT:** Please fix the guide cover of the threaded pin securely during this procedure in order to prevent the integrated pressure spring from unintentional extension.

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In step 4, please remove the "guide cover" [A] with care. In doing so, fix the "tensioning element" [B]. Remove the "pressure spring" [C] first and the "tensioning element" [B] subsequently. Also remove the "O-ring gaskets" [D] from the holders. During the assembly these have to be inserted prior to the cross pins in order to guarantee a precise holding.

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

The components can now be cleaned, e.g. in the **wash basket REF 85.194.00**, as shown in the picture. The reprocessing guidance “0703091007\_WIEDERAUFBEREITUNG\_HZ\_BA.pdf” can be found on our website <https://www.zepf-dental.com/en/user-information/> and must be followed accordingly.

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Carry out assembly of the Benex® Extractor in reverse order of disassembly. All moving parts like e.g. the threaded pin or the knurled screw should be oiled with medical white oil. For this purpose we recommend the **Medical White Oil** from the company Helmut Zepf Medizintechnik GmbH.

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The Benex® Extractor may only be **sterilized in assembled condition** (German RKI guidelines).

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### RECOMMENDED CARE PRODUCTS

Medical White Oil in a 50 ml dispensing bottle REF 24.950.10 and in an Oil Pen REF 24.950.12, 18 ml

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