HELMUT ZEPF

www.zepf-dental.com

Close to you.

MANUFACTURER SINCE 1921

[www.zepf-dental.com/
anwenderinformationen](http://www.zepf-dental.com/anwenderinformationen)

CE 0483

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3-8

MD ZEPF Ratschen

Klasse: Ir/Im / UMDNS: 13-517

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CE 0483

Dear customer,

You have decided for a high quality product from our company – many thanks. In order to ensure your long-term satisfaction, we would like to explain the use of this instrument as follows:

MANUFACTURER

HELMUT ZEPF Medizintechnik GmbH
Obere Hauptstr. 16-22
78606 Seitingen-Oberflacht, Germany
Phone: +49 (0) 7464 / 9888 0
Fax: +49 (0) 7464 / 9888 88
e-mail: info@zepf-dental.com
Internet: www.zepf-dental.com

Caution

Please read the information in this leaflet thoroughly. Improper handling and care, as well as misuse, can lead to premature wear or risks for patients and users.

INTENDED PURPOSE

The torque wrench is mainly used for final insertion of dental implants and for checking sufficient primary stability. The torque wrench has an integrated ratchet function and a calibrated scale. The required torque can be set on the scale, which helps the surgeon not to tighten the implants or components too much or too little. The torque wrench is made of high-alloy stainless steel and has a core housing to accommodate various bits / adapters.

Caution

Instruments for inserting implants and especially implants may only be used by persons who are specially trained or instructed for this purpose. The implant manufacturer's specifications also have to be considered and adhered to under all circumstances.

CONTRAINDICATION

Insufficient quality of bone mass for insertion of dental implants. Infections or ulcers at the point of insertion in the bone.

UNDESIRABLE SIDE EFFECTS:

In the case of a patient's metal intolerance, contact with the instrument, which is manufactured from medical stainless steel, can still trigger hypersensitivity (an allergic reaction).

Given such a reaction, the procedure has to be discontinued and the necessary steps initiated.

Caution

These instruments may only be used by persons who are specially trained or instructed for this purpose.

APPLICATION AND SAFETY INFORMATION

Prior to initial use and before each further use, all instruments have to be completely cleaned, disinfected and sterilized!

It is very important to inspect every product prior to use for damage and wear, e.g. cracks, breaks or defects. Never use damaged products.

The torque ratchets from Helmut Zepf Medizintechnik GmbH are delivered non-sterile. The products must first be cleaned and sterilized. Compatibility with other products has to be investigated in advance.

RECONDITIONING

The ratchet was developed for sterilization in the autoclave and has proven its worth for this kind of sterilization. However, it firstly always has to be cleaned and sterilized. The ratchet is classified by the RKI (Robert Koch Institute) as "Critical B". Please observe the instructions provided under the item "Cleaning and disinfection". Compatibility with other products like implants and their abutments has to be checked in advance.

Based on the product design, the materials used and the intended purpose, no limit on the maximum number of executable reconditioning cycles can be defined.

Caution

The products must be reconditioned separately from similar or the same instruments to avoid mixing up the individual components.

MEASUREMENT ACCURACY

The measurement accuracy of the torque ratchets is ensured by the spring used during calibration.

Calibration of the ratchet was performed on the basis of the following values:

Characteristic:

Setting at 10 Ncm / tolerance: -1.1 / +1.2 Ncm
Setting at 20 Ncm / tolerance: -1.4 / +1.6 Ncm
Setting at 30 Ncm / tolerance: -1.5 / +1.8 Ncm
Setting at 40 Ncm / tolerance: -1.7 / +2.0 Ncm

The implant manufacturer's specifications always have to be observed.

GENERAL PRINCIPLES

All products must be cleaned, disinfected and sterilized before each use; this applies in particular to the first use after delivery, as all products are delivered non-sterile (cleaning and disinfection after removal of the transport protective packaging; sterilization after packaging). Thorough cleaning and disinfection is an indispensable requirement for effective sterilization.

Please note, as part of your responsibility for the sterility of the products during use that

- only suitable equipment and product-specific validated procedures are ever used for cleaning/disinfection and sterilization
- the equipment used (washer-disinfector or WD, sterilizer, etc.) are to be regularly maintained, inspected and calibrated and
- the validated parameters are observed for each cycle.

Please also adhere to the legal requirements applicable in your country as well as the hygiene provisions of the medical practice or hospital. This applies in particular to the various specifications (e.g. in Germany according to Annex 7 of the recommendation for reprocessing from KRINKO (German Commission for Hospital Hygiene and Infection Prevention), RKI (Robert Koch Institute), BfArM (German Federal Institute for Drugs and Medical Devices) regarding effective prion inactivation (not applicable for USA).

Remark:

Reconditioning may only be carried out by trained staff in the central sterilization service department of the hospital or in the reconditioning room of the medical practice. Hospitals or medical practices are also responsible for the selection and application of the necessary protective equipment and hygiene measures.

MD ZEPF Ratchets
CE 0483
PREPARATION AND TRANSPORT

The products must be stored and transported to the reprocessing site in a closed container to avoid damage to the instruments and contamination of the environment.

CLEANING AND DISINFECTION
Principles

For cleaning and disinfection, if possible an automated procedure (washer-disinfector or WD) should be used. A manual procedure - even using an ultrasonic bath - should only be used if an automated procedure is not available or in accordance with country-specific requirements (e.g. in Germany the automated procedure is mandatory for critical B-products) due to the significantly lower effectiveness and reproducibility.

Pretreatment

Coarse impurities must be removed from the products immediately after use (within 2 h maximum). If, due to the duration of use or as a result of organizational aspects, this time cannot be adhered to, the user must determine and validate measures autonomously to prevent the contamination from drying on:

1. Disassemble the products in accordance with the specific dismantling/assembly instructions.
2. Rinse the products at least 1 min. under running water (temperature < 35°C/95°F). Moving parts should be shaken back and forth at least five times during pre-rinsing. Where applicable, brush and rinse all lumens of the products at least five times.
3. Place the disassembled products in a sufficiently large pre-cleaning bath (in an ultrasonic bath that has not been activated yet) for the specified soaking time so the products are completely submerged. Ensure that the products do not touch each other. Support pre-cleaning by completely brushing all internal and external surfaces at the beginning of the soaking time.

Caution

Be careful with products with narrow gaps in which brush bristles can get stuck! Brushes for the channels must be slightly larger than the respective channel inside diameter; the shaft length of the brush must be at least as long as the channel. Moving parts should be shaken back and forth at least five times during pre-cleaning.

Where applicable: Rinse all lumens of the products at least five times at the beginning and end of the soaking time (aids and minimum volume depending on the cavity to be rinsed).

4. Activate the ultrasound for another minimum soaking time (but not less than 5 min.).
5. Then remove the products from the pre-cleaning bath and rinse thoroughly (at least 1 min.) with water at least five times. Moving parts should be shaken back and forth at least five times during post-rinsing. Where applicable: Rinse all lumens of the products at least five times (aids and minimum volume depend on the cavity to be rinsed).

When selecting the cleaning agent² to be used, please take into consideration that

- they are generally suitable for cleaning invasive medical devices made of metals and plastics
- the cleaning agent is suitable for ultrasonic cleaning (no foam formation)
- the cleaning agent is compatible with the products (see chapter "Material Stability").

The concentrations, temperatures and soaking times, as well as rinsing instructions specified by the manufacturer of the cleaning agent or cleaning and disinfecting agent, must be strictly adhered to. Only use freshly prepared solutions, only sterile or low-bacteria (max. 10 microbes/ml) and low endotoxin (max. 0.25 endotoxin units/ml) water (e.g. purified water/highly purified water¹) and only a soft, clean and lint-free cloth and/or filtered air for drying.

Caution

Take care in case of products with rough surfaces, threads, sharp edges or comparable aspects on which particles from the cloth can be caught!

MECHANICAL CLEANING / DISINFECTION (washer-disinfector or WD)

When selecting the WD, ensure that

- the WD generally has verified effectiveness (e.g. The German Society for Hygiene and Microbiology (DGHM) or FDA approval/cleanance/registration or CE marking according to DIN EN ISO 15883)
- if possible, a tested program for thermal disinfection (A0 value > 3000 or - for older devices - at least 5 min. at 90°C/194°F) is used (in chemical disinfection, danger of disinfecting agent residues on the products)
- the program used is suitable for the products and includes sufficient rinsing steps (at least three degrading steps after cleaning (or neutralisation if applied) or conductance-based rinsing control recommended to effectively prevent detergent residues)
- only sterile or low-bacteria (max. 10 bacteria/ml) and low endotoxin (max. 0.25 endotoxin units/ml) water (e.g. purified water/highly purified water¹) is used
- the air used for drying is filtered (oil-free, low-bacteria and low-particle) and
- the WD is regularly maintained, checked and calibrated.

When selecting the cleaning system to be used, ensure that

- it is generally suitable for cleaning invasive medical devices made of metals
- providing no thermal disinfection is used - a suitable disinfecting agent with verified effectiveness (e.g. The German Association for Applied Hygiene (VAH) / The German Society for Hygiene and Microbiology (DGHM) or FDA/EPA approval/cleanance/registration or CE marking) is also used and that this is compatible with the cleaning agent used and
- the chemicals used are compatible with the products (see chapter "Material Stability").

The concentrations, temperatures and soaking times, as well as rinsing instructions specified by the manufacturer of the disinfecting agent, must be strictly adhered to.

Procedure:

1. Place the disassembled products coming from the pre-cleaning into the WD. Ensure that the products do not touch each other.
2. Start the program.
3. Remove the products from the WD at the end of the program
4. Inspect and pack the products as soon as possible after removal (see chapter "Inspection" and "Packaging", possibly after additional drying in a clean area).

Verification of the products' general suitability for effective automated cleaning and disinfecting was provided by an independent, officially accredited and recognized (§ 15 (5) German Medical Devices Act (MPG)) test laboratory using the WD G 7836 CD (thermal disinfection, Miele & Cie. KG, Gütersloh) and the neodisher MediClean forte pre-cleaning and cleaning agent (Dr. Weigert GmbH & Co. KG, Hamburg). Here, the procedure described above was taken into consideration.

MD ZEPF Ratchets
CE 0483
MANUAL CLEANING AND DISINFECTION

When selecting the cleaning and disinfecting agent to be used, ensure that

- it is generally suitable for cleaning and disinfecting invasive medical devices made of metals
- the cleaning agent is suitable for ultrasonic cleaning (no foam formation)
- a suitable disinfecting agent with verified effectiveness (e.g. The German Association for Applied Hygiene (VAH) / The German Society for Hygiene and Microbiology (DGHM) or FDA/EPA approval/cleanance/registration or CE marking) is also used and that this is compatible with the cleaning agent used and
- that the chemicals used are compatible with the products (see chapter "Material Stability").

Combined cleaning/disinfecting agents should not be used if possible.

The concentrations, temperatures and soaking times specified by the manufacturer of the cleaning and disinfecting agent, must be strictly adhered to. Only use freshly prepared solutions, only sterile or low-bacteria (max. 10 bacteria/ml) and low endotoxin (max. 0.25 endotoxin units/ml) water (e.g. purified water/highly purified water)¹ and only a soft, clean and lint-free cloth and/or filtered air for drying.

Caution

Be careful with products with rough surfaces, sharp edges or comparable aspects on which particles from the cloth can be caught!

Procedure: Cleaning

1. Place the disassembled products coming from pre-cleaning in a sufficiently large cleaning bath (in an ultrasonic bath that has not been activated yet) for the specified soaking time so that the products are completely submerged. Ensure that the products do not touch each other. Support the cleaning process by completely brushing all internal and external surfaces with a soft plastic brush and shake moving parts back and forth several times.
2. Activate the ultrasound for another minimum soaking time (but not less than 5 min.).
3. Then remove the products from the cleaning bath and rinse thoroughly (for at least 1 min.) with water at least five times.
4. Inspect the products (see chapter "Inspection").

Disinfection

5. Place the products for the predefined soaking time in the cleaning bath so that the products are completely submerged. Ensure that the products do not touch each other.
6. Then remove the products from the disinfection bath and rinse thoroughly (for at least 1 min.) with water at least five times.
7. Dry the products by blowing off/out with filtered compressed air.
8. Pack the products as soon as possible after removal (see chapter "Packaging", possibly after additional drying in a clean area).

Inspection

Check all products after cleaning or cleaning/disinfection for corrosion, damaged surfaces, chippings, contaminants and stains and discard damaged products. If the scale on the adjusting nut is no longer legible and/or it is no longer possible to adjust the adjusting nut continuously, the product must also be discarded. Products that are still contaminated must be cleaned again and disinfected.

Maintenance

Reassemble disassembled products (see specific dismantling/assembly instructions).

In the case of oiling joints, care must be taken to ensure that only instrument oils (white oil, without further additives) are used which - taking into consideration the maximum sterilization temperature applied - are approved for steam sterilization and have tested biocompatibility, and that only a small amount is applied to the joints.

Packaging

If necessary, sort the cleaned and disinfected products into the corresponding sterilization tray.

Please pack the products or the sterilization trays in sterilization containers or very large products in single-use sterilization packaging (single or double packaging) that meet the following requirements (material/process):

- DIN EN ISO/ANSI AAMI ISO 11607 (for USA: FDA clearance)
- suitable for steam sterilization (temperature stability up to at least 138°C (280°F) sufficient steam permeability)
- sufficient to protect the products or sterilization packaging from mechanical damage
- regular maintenance according to the manufacturer's specifications (sterilization container)
- do not exceed a maximum weight of 10 kg per pack / content of the sterilization container.

STERILIZATION

For sterilization, only the following sterilization methods may be used; other sterilization methods are not allowed.

Steam sterilization

- Fractionated vacuum process^{3,4} (with sufficient product drying⁵)
- Steam sterilizer in accordance with DIN EN 13060/DIN EN 285 or ANSI AAMI ST79 (for USA: FDA clearance)
- Validation in accordance with DIN EN ISO 17665 (valid IQ/OQ (commissioning) and product-specific performance assessment (PQ))
- Maximum sterilization temperature 134°C (273°F); plus tolerance in accordance with DIN EN ISO/ANSI AAMI ISO 17665
- Sterilization time (soaking time at sterilization temperature):

Country	fractional vacuum procedure
Germany	at least 5 min. ⁶ at 134°C (273°F)
USA	min. 4 min. at 132°C (270°F), drying time min. 20 min ⁴
France	at least 5 min. ⁶ at 134°C (273°F) if required for prion inactivation sterilization time 18 min.
other countries	at least 5 min. ⁶ at 132°C (270°F) / 134°C (273°F)

Gravity displacement is not recommended⁴.

Verification of the general suitability of the products for effective steam sterilization has been provided by an independent, officially accredited and recognized (§ 15 (5) MPG) test laboratory using the HST 6x6x6 steam sterilizer (Zirbus technology GmbH, Bad Grund) and using the fractionated vacuum procedure and LAWTON ME-DOIL (oils for joints and friction surfaces). Here, typical conditions in clinics and medical practices and the procedure described above were taken into consideration.

The flash sterilization process is generally not permitted.

Do not use hot air sterilization, radiation sterilization, formaldehyde or ethylene oxide sterilization, or plasma sterilization.

MD ZEPF Ratchets
CE 0483
Material Stability

When selecting cleaning and disinfecting agents, please ensure that they do not contain the following components:

- organic, mineral and oxidizing acids (minimum permissible pH value 5.5)
- Alkalis/strong alkaline solutions (neutral/enzymatic (max. permissible pH 8.5, mandatory requirement for products made of aluminum or other alkali-sensitive materials, see chapter "Special Instructions") or alkaline cleaner (max. permitted pH 11, mandatory requirement for products with intended use in prion-critical areas, e.g. In accordance with Appendix 7 of the KRINKO RKI BfArM recommendation)
- Organic solvents (e.g. alcohols, ethers, ketones, benzines)
- Oxidizing agents (e.g. hydrogen peroxide)
- Halogens (chlorine, iodine, bromine)
- Aromatic/halogenated hydrocarbons

Never clean products, sterilization trays or sterilization containers with metal brushes or steel wool.


All products, sterilization trays and sterilization containers may only be exposed to temperatures under 138°C (280°F).

STORAGE

After sterilization, the products must be stored dry in the sterilization packaging dust-free at room temperature. Avoid large temperature fluctuations. This prevents the formation of condensate and the related corrosion. In the case of proper storage, there is no shelf-life specification or functional restriction for instruments supplied non-sterile after production.

Caution

When storing the instrument, the adjusting screw must always be screwed in up to the first mark. Permanent pressure on the spring can change the accuracy of the torque ratchet.

REPAIRS

Repairs may only be undertaken by the manufacturer, otherwise every warranty for function expires.

BRIEF DESCRIPTION

The ratchet head is lockable to disable the torque limitation. In the case of non-lockable ratchet heads, the torque can be set with the setting screw. The torque is infinitely adjustable from 10 Ncm to 40 Ncm. With the aid of various inserts, the torque ratchet is used with different inserts to fix the implants in their final position, to screw in the locking screws, gingiva formers and titanium abutment screw.

APPLICATION

The torque ratchet is used with various inserts depending on the application, either with or without torque limitation. Prior to initial use and after every use on the patient, the torque ratchet and its inserts has to be disassembled, cleaned, dried and sterilized after assembly. See "Cleaning and Disinfection" and "Sterilization of the Instruments" in these Instructions for Use.

Important:


After each use, the setting screw has to be unscrewed to the lowest marking to relieve the spring!

We recommend annual calibration of the torque ratchet, which can be performed by Helmut Zepf Medizintechnik GmbH.

The user is responsible for performing the annual calibration at the manufacturer's site or an accredited test laboratory with suitable testing equipment.

USE WITHOUT TORQUE LIMITATION

Adjust the setting screw to 10 Ncm with several turns.

Pull the ratchet head out of the ratchet handle, turn 90° and click into place on the ratchet handle. The ratchet head can now no longer be bent.

After attaching the appropriate inserts (accessory components for manual surgical instruments) the torque ratchet can be used without the torque function.

MD ZEPF Ratchets
CE 0483
USE WITH TORQUE LIMITATION

Set the ratchet head such that it bends under lateral load. Once the ratchet head is locked, the setting screw has to be set to 10 Ncm with a few turns. Then the ratchet head is pulled out of the ratchet handle, turned 90° and clicked into place on the ratchet handle.

Tighten the setting screw until the required torque is visible on the measurement line and coincides with the marking on the ratchet handle. Note: 6 turns on the setting screw corresponds to a difference of 10 Ncm.

Abutments, accessories, etc. can now be screwed in according to the implant manufacturer's specifications until the ratchet head bends. This is a sign for reaching the set torque. Do not continue screwing in the bent state.

DISMANTLING/ASSEMBLY INSTRUCTIONS
Disassembly


Remove inserts from the ratchet head. Completely screw out the setting screw from the ratchet head.

Pull off the ratchet head from the ratchet handle and remove the torque spring.

Assembly

Push the torque spring into the ratchet handle.

Oil the contact and friction surfaces with medical white oil art. no. 24.950.10.

Push the ratchet head with the threaded rod into the ratchet handle.

Oil the joints and threads with medical white oil art. no. 24.950.10. Wipe off excess oil with a lint-free cloth.

Screw in the setting screw onto the ratchet handle down to the lowest marking.

HANDLING

Always treat all instruments with the greatest care while transporting, cleaning, care, sterilization and storage. This particularly applies for sensitive areas.

GUARANTEE

HELMUT ZEPF Medizintechnik GmbH supplies only tested and defect-free products to its customers. All of our products are designed and manufactured to comply with the most stringent quality requirements. We assume no liability for products that have been altered compared to the original, used for purposes for which they were not intended or used improperly. Products that are repaired or altered at any company other than HELMUT ZEPF or a company commissioned for this purpose by HELMUT ZEPF are no longer authorized to bear the CE symbol.

Note regarding return shipments: Due to legal regulations and, above all, to protect our employees, a signed proof of decontamination is required for each return shipment. This form must be completed, signed and attached to all product returns (complaints / repairs / other reasons for return).

Please pack the products in a way to avoid any risk of injury for our staff when unpacking.

The form is available for download on our website „https://www.zepf-dental.com/IFU/080300008_GB.pdf“.

Caution

All serious incidents related to the device must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

RECOMMENDED CARE PRODUCTS

Medical white oil in the 50 ml dosage bottle and in the oil pen, 12 ml

Physiologically harmless according to the German pharmacopoeia (DAB 10).

24.950.10

24.950.12
RELEASE VERSION - SEE FOOTER

¹ Should you consider a lower water quality as sufficient against the background of national recommendations (e.g. in Germany KRINKO/RKI/BfArM recommendation for treatment), this is your sole responsibility.

² If you use a cleaning and disinfecting agent for this purpose - e.g. for occupational safety reasons - please note that it should be aldehyde-free (otherwise fixation of blood contamination), have a tested effectiveness (e.g. VAH/DGHM or FDA/EPA approval/cleanance/registration or CE marking), be suitable for disinfecting the products and be compatible with the products (see chapter "Material Stability"). Please note that the disinfecting agent used during pre-treatment is for personal protection only and cannot replace the later disinfection step which must be carried out after cleaning.




³ at least three vacuum steps

⁴ The use of the less effective gravity displacement is only permissible if the fractionated vacuum procedure is not available, it requires significantly longer sterilization times and must be validated by the user for specific products, devices, processes and parameters.

⁵ The actual required drying time depends directly on parameters that are the sole responsibility of the user (loading configuration and density, sterilizer status, ...) and must, therefore, be determined by the user. Nevertheless, drying times should not be under 20 min.

⁶ or extended sterilization time (e.g. 18 min.) for prion inactivation according to national regulations (not relevant for USA)

Explanation of symbols

	Follow electronic instructions for use		Unique device identifier
	Manufacturer		Quantity
	Date of manufacture		HIBC Barcode
	Catalog number		CE marking
	LOT no. (batch no.)		Caution
	Non-sterile		Medical device

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