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**MD Titanium pins****REF 47.560.03 - 47.561.06 Class: IIb R8 / UMDNS: 16-085****and accessories, instruments for insertion / removal****REF 47.520.01 / 47.520.07 / 47.847.12**

CE 0297

**Dear customer,**

You have decided for a high-quality product from our company – many thanks.

In order to ensure your long-term satisfaction, we would like to explain the use of these instruments as follows:

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Internet: www.zepf-dental.com

**Attention**

Please read the information in this user information thoroughly. Improper handling and care, as well as misuse, can lead to premature wear or risks for patients and users. Please also note the inscriptions

on the label.

**INTENDED PURPOSE / INDICATION**

The intended purpose is the osteosynthetic restoration of ridge defects in dentistry. The titanium pin is hammered into the bone to fix membranes or titanium meshes, to isolate the bone from gingival soft tissue during regeneration and/or to fix the augmented bone / bone graft substitute which consists of autologous/ allogeneic bone material, or a mixture.

For insertion we recommend a pin applicator **REF** 47.520.01 or 47.520.07 and for removal a pin remover **REF** 47.847.12 from Helmut Zepf Medizintechnik GmbH.

**Contraindication**

An insufficient amount or quality of bone mass to allow stable fixation. Infection or ulcers where the titanium pin is to be placed. Material hypersensitivity, i.e. patient reactions to foreign bodies.

**Attention**

Instruments and especially titanium pins may only be used by persons who are specially trained or instructed for this purpose.

It is recommended to explant the titanium pins after the recovery of the alveolar ridge defects using a pin remover art. no. 47.847.12 from Helmut Zepf Medizintechnik GmbH. Upon evaluation by the physician, the physician can decide under his sole responsibility whether the titanium pins can remain in the patient's body. Care must be taken to ensure that the Titanium Pins are not swallowed or aspirated by the patient.

**APPLICATION & SAFETY INFORMATION**

It is very important to inspect every product prior to use for visible damage and wear, e.g. cracks, breaks or defects. Never use damaged products.

**SURFACE OF THE TITANIUM PINS**

The surface of all titanium pins is chemically passive; they are not magnetic.

**CORRECT SELECTION OF THE TITANIUM PINS**

We refer to the fact that titanium pins can only fulfill their function if the following basic rules are observed:

- Severe deformations of the titanium pins have to be avoided.
-  The reuse of titanium pins is not permitted.
- It is absolutely necessary to inform the patient of the advantages and disadvantages of titanium pins!

**RECONDITIONING**

The titanium pins themselves may only be used once under any circumstances. As soon as an application is concluded, the titanium pins must not be used again.

Following treatment of a patient infected with Creutzfeldt Jakob disease, the instruments and accessories for insertion / use must not be reused or reconditioned either. The instruments have to be disposed of.

**GENERAL PRINCIPLES**

All products must be cleaned, disinfected and sterilized before each use; this particularly applies to the first use after delivery, as all products are delivered unsterilized (cleaning and disinfection after removing the transport protection packaging; sterilization after packaging). A thorough cleaning and disinfection is an indispensable requirement for effective sterilization.

Please note as part of your responsibility for the sterility of the products during use that

- generally, only suitable equipment and product-specific validated procedures are to be used for cleaning/disinfection and sterilization,
- the equipment used (WD, sterilizer, etc.) are to be regularly maintained and inspected, and
- the validated parameters are to be observed for each cycle.

Please also adhere to the legal requirements applicable in your country as well as the hygienic requirements of the medical practice or hospital. This particularly applies to the various requirements (e.g. in Germany according to attachment 7 of KRINKO RKI BfArM recommendation for processing) regarding to an effective prion inactivation (not applicable for the USA).

**Remark:**

Processing must be performed only by qualified staff in the central sterilization service department of the hospital or in the processing room of the medical practice. Hospital or medical practice are responsible for selection and application of required protective equipment and hygienic measures.

**PREPARATION AND TRANSPORT**

The titanium pins must be stored and transported to the preparation site in a closed container to avoid damage to the titanium pins and/or the packaging and contamination of the environment.

**MD Titanium pins**

CE 0297

**CLEANING AND DISINFECTION****Principles**

For cleaning and disinfection, if possible an automated procedure [WD (washer-disinfector)] should be used. A manual procedure – even using an ultrasound bath – should only be used according to country specific requirements (e.g. in Germany for critical B products automated procedure binding) and if an automated procedure is not available due to the significantly lower effectiveness and reproducibility.

**AUTOMATED CLEANING AND DISINFECTION [WD (CLEANING AND DISINFECTION DEVICE)]**

When selecting the WD, ensure that

- the WD generally has verified effectiveness (e.g. DGHM or FDA approval/clearance/registration or CE marking in accordance with DIN EN ISO 15883),
- if possible, a tested program for thermal disinfection (A0 value  $\geq 3000$  or – for older devices – at least 5 min at 90 °C/194 °F) is used (in chemical disinfection danger of disinfecting agent residues on the products),
- the program used is suitable for the products and contains sufficient rinsing steps (at least three degrading steps after cleaning (respectively neutralization, if applied) or conductance based rinsing control recommended in order to prevent effectively remnants of the detergents),
- for rinsing only sterile (max. 10 bacteria/ml) and low-endotoxin (max. 0.25 endotoxin units/ml) water (e.g. purified water/highly purified water)<sup>1</sup> is used,
- air used for drying is filtered (oil-free, low-bacteria and low-particle) and
- the WD is regularly maintained, inspected, and calibrated.

When selecting the cleaning system, ensure that

- it is generally suitable for cleaning medical instruments made of metals,
- providing no thermal disinfection is used – a suitable disinfecting agent with verified effectiveness (e.g. VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking) is also used and that it is compatible with the cleaning agent used, and
- the chemicals used are compatible with the products (see “Material Stability” section).

The concentrations, temperatures and soaking times specified by the manufacturer of the cleaning agent and, if applicable, the disinfecting agent as well as specifications for rinsing must be adhered to.

**Procedure:**

1. Place the products into the WD. Ensure that the products do not touch.
2. Start the program.
3. Remove the products after the program has completed.
4. Inspect and pack the products as soon as possible after removal (see “Inspection” and “Packaging” chapters, possibly after additional drying in a clean area).

The verification of products’ general suitability for effective automated cleaning and disinfecting was provided by an independent, governmentally accredited and recognized (§ 15 (5) MPG) test laboratory using the G 7836 CD washer-disinfector (thermal disinfection, Miele & Cie. GmbH & Co., Gütersloh) and the Neodisher MediClean forte pre-cleaning and cleaning agent (Dr. Weigert GmbH & Co. KG, Hamburg). Here, the procedure described above was taken into consideration.

**MANUAL CLEANING AND DISINFECTION**

When selecting the cleaning and disinfecting agent, ensure that

- it is generally suitable for cleaning and disinfecting medical instruments made of metals,
- the cleaning agent is suitable for ultrasound cleaning (no foam formation),
- a disinfecting agent with verified effectiveness (e.g. VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking) is used and that this is compatible with the cleaning agent used, and
- the chemicals used are compatible with the products (see “Material resistance” chapter).

Combined cleaning/disinfecting agents should not be used if possible.

The concentrations, temperatures and soaking times specified by the manufacturer of the cleaning and disinfecting agent as well as specifications for rinsing, must be adhered to. Only use freshly prepared solutions, sterile or low-germ (max. 10 bacteria/ml) and low-endotoxin (max. 0.25 endotoxin units/ml) water (e.g. purified water/highly purified water)<sup>1</sup> or only with a soft, clean, lint-free cloth (Attention: caution in case of products with rough surfaces, sharp edges or comparable aspects with danger of attachment of particles from the cloth!) and/or filtered air to dry.

**Procedure: Cleaning**

1. Place the products for the predefined soaking time in the cleaning bath (in an ultrasound bath that is not already activated), so that the products are completely submerged. Ensure that the products do not touch. Support the cleaning by completely brushing all internal and external surfaces with a soft plastic brush.
2. Activate the ultrasound for an additional minimum exposure time (but not less than 5 min).
3. Then remove the products from the pre-cleaning bath and rinse them at least three times thoroughly (for at least 1 minute) with water.
4. Inspect the products (see “Inspection” chapters).

**Disinfection**

5. Place the inspected products in the disinfection bath for the predefined soaking time so that the products are completely submerged. Ensure that the products do not touch.
6. Then remove the products from the disinfection bath and rinse them at least five times thoroughly (for at least 1 minute) with water.
7. Dry the products with filtered compressed air.
8. Pack the products as soon as possible after removal (see “Packaging” section, possibly after additional drying in a clean area).

The proof of the general suitability of the products for effective manual cleaning and disinfecting was provided by an independent, governmentally accredited and respected (§ 15 (5) German Law for Medical Devices) test laboratory using the Cidezyme/Enzol pre-cleaning and cleaning agent and the Cidex OPA disinfecting agent (Johnson & Johnson GmbH, Norderstedt). Here, the procedure described above was taken into consideration.

**Inspection**

Check all products after cleaning or cleaning/disinfecting for corrosion, damaged surfaces, chippings, contaminants and stains as well as remove damaged products. Any products that are still contaminated must be cleaned again and disinfected.

## MD Titanium pins

### Packaging

Sort the cleaned and disinfected products into the corresponding sterilization tray.

Please pack the products or the sterilization trays in sterilization containers in accordance with the following requirements (material/process):

- DIN EN ISO/ANSI AAMI ISO 11607 (for USA: FDA clearance)
- suitable for steam sterilization (temperature stability up to at least 138 °C (280 °F) sufficient steam permeability)
- sufficient to protect the products or sterilization packaging from mechanical damage
- undergo regular maintenance according to the manufacturer's specifications (sterilization containers)
- do not exceed a maximum weight of 10 kg per package/contents of the sterilization container.

### STERILIZATION

For sterilization, only the following sterilization methods may be used; other sterilization methods are not allowed.

#### Steam sterilization

- Fractionated vacuum procedure<sup>2,3</sup> (with sufficient product drying<sup>4</sup>)
- Steam sterilizer in accordance with DIN EN 13060/DIN EN 285 or ANSI AAMI ST79 (for USA: FDA clearance)
- Validated in accordance with DIN EN ISO 17665 (valid IQ/OQ (commissioning) and product-specific performance assessment (PQ))
- maximum sterilization temperature 134 °C ((273 °F); plus tolerance in accordance with DIN EN ISO 17665)
- Sterilization time (exposure time at sterilization temperature):

Country	fractionated vacuum procedure
Germany	at least 5 min <sup>5</sup> at 134 °C (273 °F)
USA	at least 4 min at 132 °C (270 °F), drying time at least 20 min <sup>3</sup>
France	at least 5 min <sup>5</sup> at 134 °C (273 °F) if required for prion inactivation, sterilization time 18 min
other countries	at least 5 min <sup>5</sup> at 132 °C (270 °F) / 134 °C (273 °F)

Gravity displacement is not recommended<sup>3</sup>.

Verification of the general suitability of the products for effective steam sterilization was provided by an independent, governmentally accredited and respected (§ 15 (5) MPG) test laboratory using the HST 6x6x6 steam sterilizer (Zirbus technology GmbH, Bad Grund) and using the fractionated vacuum procedure, as well as the instrument oil LAWTON MEDOIL. Here, the typical conditions in the clinic and medical practice and the procedure described above were taken into consideration.

The flash sterilization procedure is generally not permitted.

Do not use dry heat sterilization, radiation sterilization, formaldehyde or ethylene oxide sterilization, or plasma sterilization.

### Material Stability

When selecting the cleaning and disinfection agents, please ensure that they do not contain the following components:

- Organic, mineral and oxidizing acids (minimum permitted pH value 5.5)
- Alkalis/strong alkalis (neutral/enzymatic (max. permitted pH 8.5, mandatory requirement for products made of aluminum or other alkali-sensitive materials, see "Special Instructions" section) or alkaline cleaner (max. permitted pH 11, mandatory requirement for products with intended application in prion-critical areas, e.g. in accordance with Annex 7 of KRINKO RKI BfArM recommendation for treatment) recommended)
- Organic solvents (e.g. alcohols, ethers, ketones, benzines)
- Oxidizing agents (e.g. hydrogen peroxides)
- Halogens (chlorine, iodine, bromine)
- Aromatic/halogenated hydrocarbons

Never clean products, sterilization trays or sterilization containers with metal brushes or steel wool.

All products, sterilization trays and sterilization containers can only be exposed to temperatures under 138 °C (280 °F).

### STORAGE

Please store the titanium pins in a dry, clean, and dust-free environment. Avoid large temperature fluctuations. This prevents the formation of condensate and related corrosion. In case of proper storage, there is no shelf-life specification or functional restriction for titanium pins supplied non-sterile after production.

After sterilization, the products must be stored dry and free of dust in the sterilization packaging.

### REPAIRS

Repairs are not permissible under any circumstances.

### INFORMATION ON VALIDATION OF RECONDITIONING

It is the responsibility of the user to validate their procedure appropriately.

### HANDLING

Always treat all instruments and titanium pins with the greatest care in transporting, cleaning, care, sterilization and storage. This particularly applies for sensitive areas.

### WARRANTY

Our company exclusively delivers inspected and flawless products to its customers. All our products are designed and manufactured such that the highest quality standards are fulfilled. Liability for products that are modified as compared with the original, misused or improperly treated or used, is excluded.

**HELMUT ZEPF MEDIZINTECHNIK GMBH ACCEPTS NO LIABILITY IF THIS CUSTOMER INFORMATION HAS BEEN VERIFIABLY VIOLATED.**

**THESE PRODUCTS ARE ONLY INTENDED FOR SINGLE USE.**

**MD Titanium pins**
**CE 0297**
**DECONTAMINATION INFORMATION**

On the basis of statutory regulations and above all to protect our employees, we need a signed "Evidence of decontamination" with every return.

Please ensure that this "Evidence of decontamination" is enclosed filled out and signed with **every product return** (complaint / repair / other reason for return) and package the product such that there is no risk of injury for our personnel in goods receiving during unpacking.

**This form is available for download on our website "www.zepf-dental.com".**

**Attention**

Any serious-incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

**A reporting form for this purpose can be found on our website.**

**RELEASE VERSION – SEE FOOTER**

<sup>1</sup> In case of consideration of a lower water quality as sufficient based on the background of national recommendations (e.g. in Germany KRINKO RKI BfArM recommendation for processing).

<sup>2</sup> at least three vacuum steps

<sup>3</sup> The use of the less effective gravity displacement is only permitted if the fractionated vacuum procedure is not available. It requires significantly longer sterilization times and must be validated by the user for each specific product, device, procedure and parameter.

<sup>4</sup> The actual required drying time depends directly on the parameters, which are the sole responsibility of the user (loading configuration and density, sterilizer status) and must, therefore, be determined by the user. Nonetheless, drying time should not be under 20 min.

<sup>5</sup> or 18 min (prion inactivation, not relevant for the USA)

**Explanation of symbols**

	Observe the electronic instructions for use		Unique device identifier
	Manufacturer		Quantity
	Date of manufacture		HIBC Barcode
	Catalog number		CE marking
	LOT no. (batch no.)		Caution
	Non-sterile		Medical device

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