Trephinen

REF 08.910.XX, 08.911.XX, 08.912.XX, 08.914.XX

Klasse: IIa R6 / UMDNS (Universal Medical Device Nomenclature System): 14-147

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Trephines

**REF 08.907.XX, 08.910.XX, 08.911.XX,
08.912.XX, 08.913.XX, 08.914.XX**

Class: IIa R6 / UMDNS (Universal Medical Device Nomenclature System): 14-147

Dear customer,

You have decided for a high quality product from our company – many thanks.
In order to ensure your long-term satisfaction, we would like to explain the use of these instruments as follows:

MANUFACTURER

HELMUT ZEPF Medizintechnik GmbH
Obere Hauptstr. 16-22
78606 Seitingen-Oberflacht, Germany
Phone: +49 (0) 7464 / 9888 0
Fax: +49 (0) 7464 / 9888 88
e-mail: info@zepf-dental.com
Internet: www.zepf-dental.com

The trephines from Helmut Zepf Medizintechnik GmbH are supplied non-sterile. The products have to be firstly cleaned and sterilized. Compatibility with other products has to be investigated in advance.

RECONDITIONING

Based on the product design, the materials used and the intended purpose, no limit of a maximum number of executable reconditioning cycles is defined.

PREPARATION AND TRANSPORT

The instruments must be stored and transported to the preparation site in a closed container to avoid damage to the instruments and contamination of the environment.

Attention

Please read the information in this leaflet thoroughly. Improper handling and care, as well as misuse, can lead to premature wear or risks for patients and users. Please also note the inscriptions on the packaging.

INTENDED PURPOSE / INDICATION

Processing and machining autologous bone internally (surgically invasive) and externally (extracted bone) outside the patient. Ring shaped / cylindrical bone blocks can be extracted with the trephines. A central hole allows these to be directly implanted at another recipient site or further processed in a bone mill. The instrument can be operated manually with a handle or used with a rotating handpiece switched on.

CONTRAINDICATION

An insufficient amount or quality of bone mass to allow extraction of bone / bone rings.
Infections or ulcers at the extraction site.

UNDESIRABLE SIDE EFFECTS

In the case of a patient's metal intolerance, contact with the instrument, which is manufactured from medical stainless steel, can still trigger a hypersensitivity (allergic) reaction. Given such a reaction, the procedure has to be discontinued and the necessary steps initiated.

Attention

These instruments may only be used by persons who are specially trained or instructed for this purpose.

APPLICATION AND SAFETY INFORMATION

All instruments must be completely cleaned, disinfected and sterilized before they are put into use for the first time and prior to each further use!

It is very important to inspect every product prior to use for damage and wear, e.g. cracks, breaks or defects. Never use damaged products.

MACHINE CLEANING

Place the instruments in a strainer bowl on the rack cart or on the inserts of the cart and start the cleaning process.

01. 1 min. pre-rinsing with cold water
02. Emptying
03. 3 min. pre-rinsing with cold water
04. Emptying
05. 5 min. washing at 55 °C (131 °F) with 0.5% alkaline cleaner or 45 °C (113 °F) with enzymatic cleaner
06. Emptying
07. 3 min. neutralization with warm tap water (< 40 °C (104 °F)) and neutralizer
08. Emptying
09. 2 min. intermediate rinsing with warm tap water (< 40 °C (104 °F))
10. Emptying

DISINFECTION

Perform machine thermal disinfection taking the national requirements concerning the A₀ value (see ISO 15883) into account.

DRYING

Drying the instruments with the drying cycle of the cleaning / disinfection device. If necessary, manual drying can also be achieved with the aid of a lint-free cloth.

Trephines

**REF 08.907.XX, 08.910.XX, 08.911.XX,
08.912.XX, 08.913.XX, 08.914.XX**

Class: IIa R6 / UMDNS (Universal Medical Device Nomenclature System): 14-147

MANUAL CLEANING

Prepare a cleaning bath according to the manufacturer's specification.

01. Rinse the product under cold tap water (< 40 °C (104 °F)) until all visible contamination has been removed. Stubborn dirt has to be removed with a soft brush.
02. Completely insert the products in the prepared cleaning bath. Comply with the exposure time in the manufacturer's specifications.
03. Manually clean the inserted instruments with a soft brush. All surfaces have to be brushed off several times.
04. Rinse the products very thoroughly under flowing tap water to remove the cleaning agent without leaving residues.

DISINFECTION

Prepare a disinfection bath according to the disinfection manufacturer's specifications.

Place the instruments in the disinfection bath and observe the suggested exposure time.

Rinse the products very thoroughly with DI water to remove the cleaning agent without leaving residues.

FUNCTIONAL INSPECTION AND PACKAGING

Visual inspection for cleanliness; possibly assembly and functional test. Repeat the reconditioning process as necessary until the instrument is visually clean.

Standard-compliant packaging of the instruments for sterilization in accordance with ISO 11607-01 and EN 868.

STERILIZATION

Sterilization of the products with the fractionated pre-vacuum process (according to DIN EN ISO 13060 / DIN EN ISO 17665) under consideration of the respective national requirements.

3 pre-vacuum phases with at least 60 mbar pressure.

Heating to a sterilization temperature of at least 132 °C (270 °F); max. 137 °C (279 °F). Shortest holding time: 3 min.

Drying time: at least 10 min.

STORAGE

Store instruments in a dry, clean and dust-free environment. Major temperature fluctuations have to be avoided. This prevents formation of condensate and the resulting corrosion. A shelf-life specification or functional restriction is not given for surgical or dental

instruments supplied non-sterile after manufacture if they are stored properly.

REPAIRS

Repairs may only be undertaken by the manufacturer, otherwise every warranty for function expires.

VALIDATION OF RECONDITIONING

It is the responsibility of the user to validate their procedure appropriately.

BRIEF DESCRIPTION OF CLEANING / CARE

Attention

Do not use any damaged instruments.

If an instrument can be disassembled, it should be cleaned in its disassembled state.

Cleaning and rinsing should be performed immediately after use, as otherwise particles or dried secretions can stick, which may impede subsequent cleaning and sterilization or render it impossible. Instruments must be thoroughly cleaned of foreign bodies. Use warm water and commercially available instrument cleaning agents. Observe the relevant instructions for use for these cleaning agents. Do not use any corrosive cleaning agents. Agents with a neutral pH value (7.0) are most suitable.

Do not use any abrasive cleaners.

Only use soft brushes.

Rinse and brush all channels and tubes.

Ultrasonic cleaning is suitable, but take care that connections loosened by the vibrations are fastened again. Suitable for disinfection in the washing machines up to 95 °C (203 °F).

Rinse carefully with distilled water.

HANDLING

Always treat all instruments with the greatest care in transporting, cleaning, care, sterilization and storage. This particularly applies for sensitive areas.

WARRANTY

Our company exclusively delivers inspected and flawless products to its customers. All our products are designed and manufactured such that the highest quality standards are fulfilled. Liability for products that are modified as compared with the original, misused or improperly treated or used, is excluded.

HELMUT ZEPF MEDIZINTECHNIK GMBH ACCEPTS NO LIABILITY IF THIS CUSTOMER INFORMATION HAS BEEN VERIFIABLY VIOLATED.

EVIDENCE OF DECONTAMINATION

On the basis of statutory regulations and, above all, to protect our employees, we need a signed "Evidence of decontamination" with every product return.

Please ensure that this "Evidence of decontamination" is filled in and enclosed with each product return (complaint / repair / other reason for return) and package the product such that no risk of injury exists for our personnel in incoming goods inspection while unpacking.

The form is available for downloading from our homepage at "www.zepf-dental.com/en-US/kontakt/return-shipment/".

RELEASE VERSION – SEE FOOTER