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REF 17.009.00 - 17.009.05

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
REF 17.009.00 - 17.009.05

Dear valued customer!

We thank you very much for having chosen a high-quality instrument from Helmut Zepf Medizintechnik GmbH.


In order to enjoy this instrument for a long time, we kindly ask you to follow our application notes.


For the best control, hold the elevator as illustrated.

Hi-Tools are surgical instruments and should be used as such.

Do not use excessive force.

5 Technique for molars


The instrument's tip is inserted into the sulcus towards the root in a most possible axial direction. Gently tap the Hi-tool with a hammer to expand the alveolar socket. Thereby, the Hi-tool is pushed forward in apical direction along the root of the tooth to be extracted.

For optimal placement of the instrument, the tips are worked out more delicate.

Note: These Hi-tools are the continuation of the successful X-tool concept. The instrument must not be used as a lever.

Rather, the Sharpey's fibers shall first be released. A gentle removal of the teeth will subsequently be achieved by **axial hammer taps**.

The lightweight Ferrozell hammer allows safe working and gives the patient a more agreeable feeling than a metal hammer.

Attention

Hi-tools dispose of very delicate working tips! Usage with leverage and non-planar pressure needs to be avoided. This can cause breakage of the instrument tips.

Attention

We cannot be held responsible for damages caused by improper use.

Cleaning and Reprocessing

For all cleaning processes, please consider our application guidance of autoclavable instruments and accessories according to DIN EN ISO 17664, in the currently valid version.

See document "0703091007_WIEDERAUFBEREITUNG_HZ_BA", downloadable from our website www.zepf-dental.com.

All medical devices indicated here bear the CE mark according to the Medical Device Directive 93/42/EEC.

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Muito obrigado por ter escolhido um instrumento de alta qualidade da marca Helmut Zepf Medizintechnik GmbH.