HELMUT ZEPF

REF 17.009.00 - 17.009.05

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REF 17.009.00 - 17.009.05


**Dear valued customer!**

We thank you very much for having chosen a high-quality instrument from Helmut Zepf Medizintechnik GmbH.

In order to enjoy this instrument for a long time, we kindly ask you to follow our application notes.

Explanation of symbols	
	Medical device
	Catalog number
	Caution
	CE marking


MADE  IN GERMANY


For the best control, hold the elevator as illustrated.

-Tools are surgical instruments and should be used as such.

Do not use excessive force.

**5 Technique for molars**


The instruments tip is inserted into the sulcus towards the root in a most possible axial direction. Gently tap the -tool with a hammer to expand the alveolar socket. Thereby, the -tool is pushed forward in apical direction along the root of the tooth to be extracted.

For optimal placement of the instrument, the tips are worked out more delicate.

**Note:** These -tools are the continuation of the successful -tool concept. The instrument must not be used as a lever.

Rather, the Sharpey's fibers shall first be released. A gentle removal of the teeth will subsequently be achieved by **axial hammer taps**.

The lightweight Ferrozell hammer allows safe working and gives the patient a more agreeable feeling than a metal hammer.

**Attention**

-tools dispose of very delicate working tips! Usage with leverage and non-planar pressure needs to be avoided. This can cause breakage of the instrument tips.

**Attention**

We cannot be held responsible for damages caused by improper use.

**Cleaning and Reprocessing**

For all cleaning processes, please consider our application guidance of autoclavable instruments and accessories according to DIN EN ISO 17664, in the currently valid version.

See document "0703091007\_WIEDERAUFBEREITUNG\_HZ\_BA", downloadable from our website [www.zepf-dental.com](http://www.zepf-dental.com).

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Muito obrigado por ter escolhido um instrumento de alta qualidade da marca Helmut Zepf Medizintechnik GmbH.