





KlingEX^{REF} 46.005.70


- 🇩🇪 **Sicherheitsklingenabnehmer**
- 🇬🇧 **Safety Blade Remover**
- 🇫🇷 **Changeur de lames de bistouri**
- 🇪🇸 **Remover hojas de bisturís seguramente**
- 🇮🇹 **Leva lama di sicurezza**



- 🇩🇪 **positionieren**
- 🇬🇧 **positioning**
- 🇫🇷 **positionner**
- 🇪🇸 **posicionar**
- 🇮🇹 **posizionare**



- 🇩🇪 **Klinge anheben durch Pressen**
- 🇬🇧 **lift blade by pressing**
- 🇫🇷 **élever la lame en pressant**
- 🇪🇸 **elevant la hoja por presionar**
- 🇮🇹 **lama levare tramite pressione**



- 🇩🇪 **abziehen**
- 🇬🇧 **pull off**
- 🇫🇷 **retirer**
- 🇪🇸 **descontar**
- 🇮🇹 **estrarre**