





newsticker

H-Tool Set

acc. to Dr. Hildebrand

The new H-Tools developed acc. to Dr. Hildebrand are a continuation of the successful X-Tool concept.

The instrument tips are very flat and sharp pointed. The instruments are inserted along the root axis and the alveolar cavity is widened by advancing the H-Tool apically along the root of the tooth to be extracted with gentle hammer blows.

The light ferrozell hammer with a diameter of 35 mm enables work to be carried out safely.

Thanks to the extremely light construction, this hammer produces a far less unpleasant feeling in the patient than using a hammer with hard plastic inserts or in the worst case with a metal head.

Content: H-Tool Set

Illustration **Article Description Order Quantity** 1 set 17.009.00 H-Tool Set acc. to Dr. Hildebrand organized in a washbasket, including a Ferrozell hammer and H-Tools, consisting of: 1 piece 41.503.01 Hammer, Ferrozell head Ø 35 mm, 100 g, total weight 180 g, 250 mm long, fitting washbasket 1 piece 17.009.01 H-Tool #1, straight, tip 2.5 mm, green ball 1 piece 17.009.02 H-Tool #2, straight, tip 3.5 mm, green ball 1 piece 17.009.05



85.195.00

85.181.04

85.181.05

1 piece



Washbasket 1/1 as rack for the remaining components

H-Tool #5, straight, tip 3.0 mm, gold colored ball



Profile, high, universal 130 mm / for 8 instruments

2 pieces

2 pieces



Profile, high, universal 81 mm / for 5 instruments

AESTHETIC IS THE RESULT



