

## Temporary Crown Remover

acc. to Dr. Hobl, in **ZEPF** *bionik* Handle



## Temporary Crown Remover

acc. to Dr. Hobl, in **ZEPF BIONIK** Handle

The temporary crown remover with **ZEPF BIONIK** handle features a short sickle shape.

The working tip can be adapted safely to the crown margin. The short working tips offer a very direct way to transmit the tensile force on the temporary crown, i.e. a slipping off is almost impossible.



**19.265.16**

Temporary Crown Remover, double-ended, in **ZEPF BIONIK** handle, turquoise-brightblue, straight for anterior teeth, curved for posterior teeth



**BIONIK**

reddot design award  
winner 2010

**BIONIK Universal Handle** made of PEEK high-tech plastic material – guarantees an ideal power transmission with formerly unknown sensitivity. The handle is available in 10 different, fresh basic colors.

The exchangeable working tips inserted in the ergonomic **BIONIK** handle offer highest economy and best tactile handling.



yellow



signal orange



red-purple



signal purple



lightred-magenta



turquoise-brightblue



cobalt-blue



yellow green



grey



black

AESTHETIC IS THE RESULT