

# Augmentation Material Applicator

in **ZEPF** *BIONIK* Handle



## Augmentation Material Applicator

in **ZEPF BIONIK** Handle

The augmentation material applicators allow an easy condensing of bone graft substitute around the implant. The working tips turned 15° oppositely allow the practitioner an easy access in all quadrants.



**19.603.10** Augmentation Material Applicator, double-ended, in **ZEPF BIONIK** handle, turquoise-brightblue, 3.2 mm wide, 1.2 mm thick / tips turned oppositely



**19.603.15** Augmentation Material Applicator, double-ended, in **ZEPF BIONIK** handle, cobalt blue, 4.5 mm wide, 1.2 mm thick / tips turned oppositely



**BIONIK**<sup>ZEPF</sup>

reddot design award  
winner 2010

**BIONIK Universal Handle** made of PEEK high-tech plastic material – guarantees an ideal power transmission with formerly unknown sensitivity. The handle is available in 10 different, fresh basic colors.

The exchangeable working tips inserted in the ergonomic **BIONIK** handle offer highest economy and best tactile handling.

AESTHETIC IS THE RESULT