

Clips for Rinsing Hose






Clips for Rinsing Hose


A major and recurring problem in all dental treatments and operations is the limited space in the patient's mouth. Quite often it is necessary to combine rinsing and suctioning. For that purpose, we created a new clip system which fits perfectly onto our aspirators.

The rinsing tube can now be fixed easily on the aspirator itself, thus allowing fast and precise working. Suctioning and rinsing, all in one step, has never been easier!




The plastic clips are re-usable and are appropriate for steam sterilization.

Clips to fix the rinsing hose on the aspirator

Illustration	Description	Order Quantity
	19.649.05 Set consisting of a guiding clip and two stabilization clips	1 set
	19.649.06 Stabilization Clip	1 piece
	19.649.07 Guiding Clip	1 piece



Fitting surgical aspirators, curved, with slot against sucking in

	19.649.15 Ø 1.5 mm	1 piece
	19.649.30 Ø 3.0 mm	1 piece
	19.649.50 Ø 5.0 mm	1 piece

AESTHETIC IS THE RESULT

