

ZEPF Augmentation Material Applicator



ZEPF Augmentation Material Applicator



The Augmentation Material Applicator supports the fractional filling of the subantral area with particulate augmentation material. To compact the augmentation material you use a sinus plugger.

The instrument with diameter of 5 mm is loaded laterally.

The augmentation material can be applied precisely with the injecting mechanism.

The instrument can be taken apart for optimal cleaning.



19.714.16

ZEPF Augmentation Material Applicator Ø 5 mm, straight, with a lateral opening to insert the augmentation material



41.868.07

Kirsch, Sinus 7, Plugger Ø 5.0 mm, Spoon 8 x 10 mm, Titanium **ZEPF**-Line Handle, double-ended, 17.5 cm



AESTHETIC IS THE RESULT