

Explorer 11/12  
delicate, in **ZEPF** **BIONIK** Handle



## Explorer 11/12

delicate, in **ZEPF BIONIK** Handle

The Explorer 11/12 with **ZEPF BIONIK** handle features delicate and pointed tips.

For subgingival root examination. To explore pockets, restorations, furcations and to diagnose proximal and cervical calculus and caries.



**24.440.00** Explorer 11/12, double-ended, in **ZEPF BIONIK** handle, signal purple, delicate, pointed tips, angles according to Gracey 11/12



**BIONIK**<sup>ZEPF</sup>

reddot design award  
winner 2010

**BIONIK Universal Handle** made of PEEK high-tech plastic material – guarantees an ideal power transmission with formerly unknown sensitivity. The handle is available in 10 different, fresh basic colors.

The exchangeable working tips inserted in the ergonomic **BIONIK** handle offer highest economy and best tactile handling.

AESTHETIC IS THE RESULT