

## ZEPF Contrast PA Probe Inserts



## Contrast PA Probe Inserts

The new **ZEPF** Contrast PA Probes are made of plastic material and dispose of a flexible working tip which adapts optimally to the anatomy of the pocket depth when measuring.

Colorstable, black markings on the white basic material guarantee a very good contrast for reading.

The sterilizable, exchangeable tips are available in different common graduations.

They are suited to determine the parodontal status and especially to be used on implants. Scratching of implant surfaces is avoided with these probes.

**The tips are reusable until they bend, the color fades or the graduation is not readable any more.**



**24.451.00** # 1

graduation 3/6/8/11 mm  
M4 x 0.5 mm, PU 12 pieces



**24.451.01** # CPG 12

graduation 3/6/9/12 mm  
M4 x 0.5 mm, PU 12 pieces



**24.451.02** # CPNG 22

graduation 2/4/6/8/10/12 mm  
M4 x 0.5 mm, PU 12 pieces




**24.451.03** # PCPG 11.5

graduation 3.5/5.5/8.5/11.5 mm  
M4 x 0.5 mm, PU 12 pieces




**24.451.06** # CNC

graduation 1 - 15 in mm steps,  
North Carolina  
M4 x 0.5 mm, PU 12 pieces



Flexible working tip



Packing Example



Assembly without a tool

**24.451.02** # CPNG 22

Periodontal Probe exchangeable, graduation:  
2/4/6/8/10/12 mm, M4 x 0.5 mm, PU 12 pieces

**26.193.15**

**BIONIK** Universal Handle single-ended

**QUICKFIX**, lightred-magenta, incl. 1 end cap  
The handle is available in 10 different colors.

AESTHETIC IS THE RESULT