

BIONIK™ Retraction-Thread Plugger



BIONIK™ Retraction-Thread Plugger

Retraction threads are placed to repress or retract the gingiva from the tooth neck.

Prior to taking an impression for the preparation of crowns, the thread is adapted around the prepared tooth and the gingiva.

The **ZEPF BIONIK** Instruments feature specially adapted shapes, helping the dentist to place the threads in a time-saving way; this is useful as threads are tending to get thinner.



BIONIK™

red dot design award
winner 2010

BIONIK Universal Handle made of PEEK high-tech plastic material – guarantees an ideal power transmission with formerly unknown sensitivity. The handle is available in 10 different, fresh basic colors.

The exchangeable working tips inserted in the ergonomic **BIONIK** handle offer highest economy and best tactile handling.



24.548.01 Universal Thread Plugger with angled working tips, micro serrated



24.548.02 Suture Applicator for distal / mesial use, micro serrated and specially adapted to anatomy



24.548.03 Suture Applicator for buccal / lingual use, micro serrated and specially adapted to anatomy



24.548.04 Universal Thread Plugger with round tips, micro serrated



AESTHETIC IS THE RESULT

