


ZEPF Diagnostic Set

Containing all necessary Instruments for diagnostics:

- Contrast Probe # CNC in **ZEPF bionik** Handle
- Probe # 8 in **ZEPF bionik** Handle
- Lipcare Mouth Mirror in **ZEPF bionik** Handle
- Ergonomic Dental Tweezer
- Washtray 1/3



24.991.40

Essential advantages:

- ergonomically perfect
- easy cleaning
- quickly replaceable working tips
- all instruments organized in one tray
- cost-efficient



24.062.24.

Lipcare Mouth Mirror / Retractor Combination
acc. to Dr. Preusse,

- optimal ergonomic working for the user
- patient-friendly, atraumatic, extensive retraction of the cheek

Content **ZEPF** Diagnostic Set



26.194.05

BIONIK Handle, single-ended with end cap and reduction insert (M4 x 0.5 to M2.5), lightred-magenta

24.449.30

8 Probe Tip, REF 24.099.08, interchangeable, with M2.5 metric thread in **BIONIK** Handle, REF 26.194.04, **QUICKFIX**, single-ended with end cap and reduction insert (M4 x 0.5 to M2.5), signal purple




24.449.35

CNC Contrast Probe, REF 24.451.06, interchangeable, graduation 1 - 15 in mm steps, North Carolina, M4 x 0.5 mm in **BIONIK** Handle, REF 26.193.12, **QUICKFIX**, single-ended with end cap, signal orange



Flexible working tip

ZEPF *care*
system



85.182.00

Washtray 1/3,
178 x 90 x 24

ZEPF Ergonomic Dental Tweezers

Our extra comfortable **HELMUT ZEPF** Ergonomic Dental Tweezers are pleasant to hold and have been specially designed for use in lengthy, otherwise tiring, applications.

The broader version and their somewhat shorter tips allow exerting pressures in excess of those of common tweezers.



22.025.03

Ergonomic Dental Tweezers, grooved tips, 15 cm

AESTHETIC IS THE RESULT