



## Sinus Instrument Set

acc. to Prof. Dr. Dr. Stiller



24.995.01

**STILLER**



**24.995.01**


Complete Set incl. 1 x washbasket with lid  
85.195.00 and partition for Sinus Instrument Set  
acc. to **Prof. Dr. Dr. Stiller**



## Sinus Instrument Set acc. to Prof. Dr. Dr. Stiller

Immediate sinus floor elevation  
with or without bone cover

The Sinus Instrument Set has been developed together with Prof. Dr. Dr. Stiller and is particularly suited for difficult maxillary sinus structure (septa, maxillary pillars, scarred mutations in change with intact maxillary sinus mucosa).




### Advantages of the system:

- Instruments which are adapted and adjusted perfectly to the anatomical conditions in the maxillary sinus.
- Flexible working ends. Instruments are pre-bended for the normal enantial anatomy.
- Two different kinds of instruments with blunt and sharp edges for dissecting mucosa on plane and rough internal surface of the maxillary sinus.




**41.822.01**

**Stiller Sinus Elevator**, blunt, 4.0 mm, bendable, Titanium **ZEPF-Line**, green, double-ended, 19.5 cm



**41.822.02**

**Stiller Sinus Elevator**, blunt, 2.5 mm, bendable, Titanium **ZEPF-Line**, green, double-ended, 19.5 cm



**41.822.03**

**Stiller Sinus Elevator**, blunt, 3.0 mm, bendable, Titanium **ZEPF-Line**, green, double-ended, 19.5 cm



**41.822.04**

**Stiller Sinus Elevator**, universal, 3.0 mm, Titanium **ZEPF-Line**, yellow, double-ended, 19.5 cm



**Surgical Aspirator**

Sinus-Line, slotted, titanium tip, acc. to Dr. Maty

**19.651.13** Ø 1.5 mm**19.651.14** Ø 3.0 mm**41.822.05****Stiller Sinus Elevator, universal, 6.0 mm, Titanium ZEPF-Line, blue, double-ended, 19.5 cm****41.822.11****Stiller Sinus Elevator, sharp, 4.0 mm, bendable, Titanium ZEPF-Line, red, double-ended, 19.5 cm****41.822.22****Stiller Sinus Elevator, sharp, 2.5 mm, bendable, Titanium ZEPF-Line, red, double-ended, 19.5 cm****41.822.33****Stiller Sinus Elevator, sharp, 3.0 mm, bendable, Titanium ZEPF-Line, red, double-ended, 19.5 cm****41.868.07****Kirsch Sinus 7, Plugger Ø 5.0 mm, Spoon 8 x 10 mm, Titanium ZEPF-Line Handle, double-ended, 17.5 cm****85.251.04**

Mixing Cup, stainless steel, with plastic lid, Ø 40 mm

**85.251.14**

Mixing Cup without plastic lid, stainless steel, Ø 40 mm

**08.906.029C**

Diamond, 029C, round, Ø 2.9 mm

**08.906.023C**


Diamond, 023C, round, Ø 2.3 mm

**08.902.031HF**

HB Drill, 031HF, round



## Sinus Instrument Set acc. to Prof. Dr. Dr. Stiller



24.995.01

## Surgical Aspirators

Curved, with slot against sucking in, 17.5 cm

*SINUS*-line

The specially bent Sinus-Line Aspirators with titanium tip were developed by Dr. Maty.



**19.651.13**

Surgical Aspirator, with slot, titanium tip, Ø 1.5 mm

**19.651.14**

Surgical Aspirator, with slot, titanium tip, Ø 3.0 mm

AESTHETIC IS THE RESULT