

Retractor "Whale Fin"
acc. to Dr. Maty



Retractor "Whale Fin"

acc. to Dr. Maty



37.448.30

Retractor "Whale Fin",
size 1 (30 mm),
sterilizable



37.448.35

Retractor "Whale Fin",
size 2 (35 mm),
sterilizable

Both sizes of the „Whale Fin“ Retractor are available as insert M4 x 0.5 or as complete instrument in a **ZEPF BIONIK** handle.

Advantages:

- good view of the operating and preparation field
- safe retraction of soft tissues and good visual control even in the lateral tooth area
- atraumatic – no punctual contact, but surface contact
- extra opening for the labial frenulum (upper jaw / lower jaw) with center cut-out in the middle of the retraction area / working surfaces
- wave-shaped working surface for a better adaptation in the upper jaw / lower jaw
- also usable in prosthetics (e.g. tooth preparation)

The ergonomic shape of the periosteal elevator and the retractor allow an easy cleaning and sterilization procedure.



BIONIK^{ZEPF}

reddot design award
winner 2010

BIONIK Universal Handle made of PEEK high-tech plastic material – guarantees an ideal power transmission with formerly unknown sensitivity. The handle is available in 10 different, fresh basic colors.

The exchangeable working tips inserted in the ergonomic **BIONIK** handle offer highest economy and best tactile handling.



37.446.06

Retractor "Whale Fin",
size 1 (30 mm),
in **ZEPF BIONIK**
Universal Handle,
PEEK, single-ended,
turquoise-brightblue



37.446.07

Retractor "Whale Fin",
size 2 (35 mm),
in **ZEPF BIONIK**
Universal Handle,
PEEK, single-ended,
yellow green

AESTHETIC IS THE RESULT