

ZEPF Retractor and Mouth Mirror
acc. to Dentist Beck



ZEPF Retractor and Mouth Mirror
acc. to Dentist Beck

The **ZEPF BIONIK** Retractor and Mouth Mirror is available in two sizes, featuring a heart-shaped, ergonomic and highly polished design.
The instrument is available in size 24 and 27 mm.

It is universally applicable as:

- Mouth mirror, cheek or tongue retractor, mucoperiosteal flap retractor.
- Its innovative heart shape allows ergonomic adaption in all quadrants.
- Due to its reflective surface, the retractor allows an optimal illumination of the surgical site.
- The large contact surfaces make the retractor very patient-friendly.



37.448.10



37.448.20

24.062.30

ZEPF BIONIK Retractor and Mouth Mirror, size 1, 24 mm,
ZEPF BIONIK Universal Handle, PEEK, single-ended, signal orange

24.062.31

ZEPF BIONIK Retractor and Mouth Mirror, size 2, 27 mm,
ZEPF BIONIK Universal Handle, PEEK, single-ended, signal purple



AESTHETIC IS THE RESULT