

Periosteal Elevator

Micro Raspatory acc. to Dr. Maty



Periosteal Elevator

Micro Raspatory acc. to Dr. Maty

- due to its slim design, the periosteal elevator can be used in all surgical interventions, particularly in **biological** dentistry, surgery and implantology
- safe exposure and retraction of the periosteum and good visual control
- lies perfectly in the hand and guarantees a safe support and guidance
- slip-resistant
- ergonomic shape for a better removal of the periosteum
- 3 working surfaces – front, right, left – allowing a fast flap opening of soft tissues / periosteum
- optimally suited for areas that are difficult to access, e.g. wisdom teeth in the upper and lower jaw, implants (distal)
- can be used as three-dimensional “curette” and periosteal elevator, e.g.
 - a) cleaning of the alveolus after tooth extraction
 - b) use in the apical region during apicoectomy
 - c) use in periodontal pockets (removal of soft tissues)
- the instrument is lightweight, has a good grip, is easy to handle and is slip-resistant
- the periosteal elevator can be fixed and placed easily with the index finger (even in areas that are difficult to access)
- claw and clover-shaped tips to serve as pusher to release the mucoperiosteal flap



41.864.80 Periosteal Elevator acc. to Dr. Maty, **ZEPF-Line**, Micro Raspatory, double-ended



41.854.80 Periosteal Elevator acc. to Dr. Maty, in **ZEPF BIONIK** Handle, yellow green, Micro Raspatory, double-ended



BIONIK

reddot design award
winner 2010

BIONIK Universal Handle made of PEEK high-tech plastic material – guarantees an ideal power transmission with formerly unknown sensitivity. The handle is available in 10 different, fresh basic colors.

The exchangeable working tips inserted in the ergonomic **BIONIK** handle offer highest economy and best tactile handling.

AESTHETIC IS THE RESULT