

## **ZEPF** ■ Drop-Control ➤ Blade Holder

Patented





Manufacturer since **1921** 



**ZEPF** ■ Drop-Control ▶

Blade Holder - Patent No. 102014101658

The ZEPF Drop-Control® Blade Holder is available in three diameters. This allows the user to choose according to his taste.

These instruments provide additional safety in daily surgical use when changing the blade. Pressing the slider triggers a mechanism that lifts the blade, pushes the blade further forward in its actuation path and thus finally, safely and in a controlled manner, ejects it into a sharps collection container. Afterwards, the ejector mechanism moves back to its initial position. This allows the blade to be changed safely, quickly and in a time-saving manner during treatment.

- Reduction of the risk of infection due to a contaminated blade thanks to safe ejection function
- Time saving when changing the blade during treatment
- No additional instrument required to remove the blade



46.007.02 Drop-control® Blade Holder Ø 8 mm Controlled dropping of the blade 46.007.08 Drop.control® Blade Holder Ø 12 mm 46:007.18 Drop-Control® Blade Holder Ø 9.5 mm

AESTHETIC IS THE RESULT





