

Tunneling Instrument with holes  
in **ZEPF bionik** Handle



## Tunneling Instrument with holes

in **ZEPF BIONIK** Handle

Tunneling Instrument with holes, for upper and lower jaw with **ZEPF BIONIK** handle. Both ends with hole Ø 0.8 mm for inserting the surgical suture, for easy insertion of the graft into the tunnel.



**46.040.20** Tunneling Instrument with holes, double-ended, in **ZEPF BIONIK** handle, grey, for upper and lower jaw, width 2.5 mm, each hole Ø 0.8 mm



**BIONIK**<sup>ZEPF</sup>

reddot design award  
winner 2010

**BIONIK Universal Handle** made of PEEK high-tech plastic material – guarantees an ideal power transmission with formerly unknown sensitivity. The handle is available in 10 different, fresh basic colors.

The exchangeable working tips inserted in the ergonomic **BIONIK** handle offer highest economy and best tactile handling.

AESTHETIC IS THE RESULT