


Bone Grafting Kit  
in **BIONIK** Handles




**BIONIK**<sup>ZEPF</sup>



reddot design award  
winner 2010

**Bone Grafting Kit**  
in **BIONIK** Handle

Ideally matched instruments in **ZEPF BIONIK** Handle for bone augmentation.  
With exchangeable inserts, thread M4 x 0.5, consisting of:



**47.848.00**

**Bone Grafting Kit**  
incl. 1/3 Washtray  
**REF** 85.182.00.



**47.848.01**  
**Graft Carrier Plugger**

Plugger/Spoon combination for contouring and inserting grafts. **BIONIK** Handle, lightred-magenta.



**47.848.02**  
**Graft Material Carrier**

Augmentation Applicator/Spoon Combination for safe insertion of grafts of different sizes at the recipient site. **BIONIK** Handle, yellow.



**47.848.03**  
**Graft and Sinus Packer**

Plan/Ball Plugger Instrument for condensing augmentation material. **BIONIK** Handle, turquoise-brightblue.



**24.532.19**  
**Condenser, titanium**

with graduation 3/6/9/12/15 mm for depth measurement and for condensing augmentation material around an implant. **BIONIK** Handle, yellow-green.

Ø 1.0 mm

Graduation:

- 3.0 mm
- 6.0 mm
- 9.0 mm
- 12.0 mm
- 15.0 mm

Ø 2.0 mm

AESTHETIC IS THE RESULT