



ZEPF Bone Block Holder with Stand
Your third hand!







ZEPF Bone Block Holder with Stand

Your third hand!

Fixed in the stand and in the chuck, the bone block can be processed precisely. The quick lock allows to release the bone block holder rapidly and easily from the stand. Subsequently, due to its ergonomic and light construction, the bone block holder is optimally qualified to fix the transplant in the surgical site.



47.956.20

ZEPF Bone Block Holder



47.956.30

ZEPF Stand for Bone Block Holder


47.956.10

Set:

ZEPF Bone Block Holder with Stand

Dismantling

Assembly in reverse order.



AESTHETIC IS THE RESULT